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al., this choice should be made in the pre-surgery stage, but the same concepts were applied by analyzing the master models and the position of the teeth. The result led to the choice of a removable prosthesis considering that the angle formed by the line passing through the incisal margin of the central incisors and the implant connection with the parallel to the occlusal plane was less than 30° (Prosthetic Profile) and the distance between the border of the anterior flange and the implant connection was about 5 mm (MucCov) (Figs. 41 and 42). The advantages of this choice are clearly the easy access for hygiene procedures, the ability to position the tooth correctly for the phonetics and esthetics required by the patient, and the ability to give suitable support to the maxillary lip.

Maxillary prosthesis design

A design was therefore set for the maxillary prosthesis with a milled bar on which to position the spherical retentions (OT Cap Micro, Rhein'83 Srl, Bologna, Italy) which with limited overall dimensions, may reduce the volumes necessary for the rehabilitation (Fig. 43). The first evaluation carried out with the masks resulting from the replicas already made it clear that their positioning would have to be on the palatal side. In fact, the vertical space available did not allow their housing on the occlusal side of the future bar. Regarding the construction of the bar, the pillars mounted in the oral cavity to support the old rehabilitations could not be removed since they were no longer in production and there was no specific removal kit available. Consequently, some Equator posts were manufactured and personalized by Rhein'83 Srl, Bologna, Italy; a bar was then manufactured on these using the Seeger system (Rhein'83 Srl, Bologna, Italy) (Figs. 44 to 48). It was possible to use this system to manage any passivation problems related to the fact that the specified Equator pillars were positioned to support the bar directly in the laboratory without taking the impression of their position.

Individualized Equator Post preparation

The Equator was manufactured using a scan of the platform of the original pillar. The attachment was positioned on this engagement taking into account the transmucosal path measured on the model and in the oral cavity. This system also allows management of the strong disparallelism present between the implant pillars.

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
42

		Mucosal Coverage (MucCov)		
		0mm	0-5mm	5mm
Prosthesis Profile (ProsthProfile)	≥45 degrees	A Fixed prosthesis Crown design		
	30-45 degrees	B Probably fixed prosthesis Hybrid design		
	≤30 degrees	C Removable overdenture or total prosthesis		