Fabbricante: RHEIN 83 S.R.L. Via E. Zago, 10 A,B,C
40128 Bologna, Italia • Tel. 051-244510 • Fax 051-245238 •
N. Verde 800-901172

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Instructions for use TITANIUM+TiN CONCAVE SPHERES

DEVICE DESCRIPTION AND INDICATIONS FOR USE: the TITANIUM+TiN CONCAVE SPHERES line is composed by nitruated titanium spheres for OT CAP NORMO, OT CAP MICRO, 2,2mm SPHERES and OT EQUATOR designed for the retentive diameter recovering of the worn-out spheres and hemi-spheres. **IMPORTANT:** these instructions for use are the most recent available; please read and store them carefully. **CONTRAINDICATIONS:** at the present state of knowledge there are no contraindications to the use of the TITANIUM+TiN CONCAVE SPHERES line; at the date of emission of these instructions for use no undesirable side effects have been reported. **CAUTION:** the choice of the attachment is demanded to the dentist or to the dental technician according to the prosthetic project. **MATERIAL:** all the materials used for the production of the TITANIUM+TiN CONCAVE SPHERES line are of quality and comply with the international reference standards. The patient has to inform the dentist about potential allergies. If allergic reactions occur search the source and remove it. In case of proven allergy to the materials used for the product production, do not use the TITANIUM+TiN CONCAVE SPHERES line. **WARNINGS AND SINGLE USE:** the TITANIUM+TiN CONCAVE SPHERES are single-use products. Their reuse compromises the mechanical functionality and affects the patient's health, such as the risk of infection through infectious material transmitted from one patient to another. RHEIN 83 S.R.L. assumes no responsibility for any damage derived from an improper use or reuse of the TITANIUM+TiN CONCAVE SPHERES line. **TRANSPORTATION, STORAGE AND EXPIRATION DATE:** do not damage the packaging during the transportation; store in a clean, dry place, in the original packaging away from sunlight and heat sources. The product does not expire. **WASHING AND STERILIZATION:** the content of the set is sold as NON-STERILE; the product must be washed in hot water with a soft brush and non-aggressive, non-corrosive cleaning solutions (avoid solutions containing chlorine). The dentist is required to rinse, dry, pack and store the prefabricated metal parts in cleared FDA 510(k) pouches (paper/plastic film) suitable for sterilizing and must sterilize them before use with the following parameters: steam sterilization with prevacuum cycles at 135 °C – 3 minutes, drying time 16 minutes, pressure 2.06 bar. It is recommended to sanitize the RETENTIVE CAPS and the TRANSPARENT INSERTER by immersing them in pure benzalkonium chloride solutions for 10 minutes, to rinse and dry them before use. The dentist is responsible of the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions to the numbers listed above or to the email: gianni.storni@rhein83.it. For more details on the use of the TITANIUM+TiN CONCAVE SPHERES line are of quality and meet the international standards of reference. For more informations about the use of TITANIUM+TiN CONCAVE SPHERES line please consult the catalog-manual and the leaflets. **MAINTENANCE AND PERIODIC CARE:** dentists have the responsibility to keep the proper functionality and retention of the attachments assuring the safety of the patient by constant maintenance. **GUIDELINES FOR THE PATIENTS:** patients are recommended to follow the indications provided by the dentist, to attend periodical controls and to perform a daily, accurate hygiene. the TITANIUM+TiN CONCAVE SPHERES line has not been evaluated for safety and compatibility in the RM environment. They have not been tested for heating, migration, or image artifacts in the RM environment. The safety of the TITANIUM+TiN CONCAVE SPHERES line in the RM environment is unknown. Scanning a patient who has this device could cause injury to the patient. **DISPOSAL:** follow the local laws for the disposal of medical devices. It is suggested to wash and sterilize the products before disposal. **CONTENT OF THE SETS:** all the components of the TITANIUM+TiN CONCAVE SPHERES line are sold both as spare parts and as part of a set.

TECHNICAL SPECS

OT CAP NORMO Ti+TiN CONCAVE SPHERES: external Ø 2.5mm, internal Ø 1.9mm.

OT CAP MICRO Ti+TiN CONCAVE SPHERES: external Ø 1.8mm, internal Ø 1.4mm.

OT EQUATOR Ti+TiN CONCAVE SPHERES: external Ø 2.5mm, internal Ø 2.1mm.

2,2mm Ti+TiN CONCAVE SPHERES: external Ø 2.2mm, internal Ø 1.55mm.

TRANSPARENT INSERTER: double sided steriliable plastic accessory available in different sizes for the four CONCAVE SPHERES.

DIAMOND STRIP HOLDER: metallic accessory with one side that allows the surface finishing of the worn-out sphere with the help of diamond strips and with the other side to be used as a cement spatula.

INSTRUCTIONS FOR USE

TRANSPARENT INSERTER: accessory made to facilitate the handling of the corresponding CONCAVE SPHERE; the two ends of the INSERTER are appropriately shaped to insert the CONCAVE SPHERE on the peduncle of the worn-out sphere (side A). **ATTENTION: insert the CONCAVE SPHERE into the INSERTER so that the diameter B of the CONCAVE SPHERE faces the worn-out sphere.** The B side of the INSERTER is used to put the CONCAVE SPHERE inside the RETENTIVE CAP present in the denture if the same is used as the positioner of the CONCAVE SPHERE.

Ti + TiN CONCAVE SPHERE: the worn-out sphere diameter should be reduced with a cylindrical diamond bur by testing, from time to time, the passive insertion of the CONCAVE SPHERE (positioned on the A side of the TRANSPARENT INSERTER) that shouldn't have excessive tolerance.

ATTENTION: do not over-reduce the worn-out sphere in order to avoid the fracture of the same or the decementation of the new CONCAVE SPHERE. Carefully degrease and dry the peduncle of the worn-out sphere and the CONCAVE SPHERE and deposit the correct amount of self-curing composite cement inside the CONCAVE SPHERE. Place the CONCAVE SPHERE on the peduncle of the worn sphere and wait for the polymerization of the cement (following the manufacturer's instructions for use) without moving the INSERTER. Once the cement set, remove the INSERTER and remove any cement excess. Place a PROTECTIVE DISK and a new RETENTIVE CAP in the corresponding METAL HOUSING on the CONCAVE SPHERE, check the dimensions inside the denture and reposition it with self-curing resin.

DIAMOND STRIP HOLDER: accessory to be used with a diamond strip inserted in the C side for the finishing of the worn-out sphere or using the D side as a spatula for the cement.

To avoid a sudden deterioration of the performances of the TITANIUM+TiN CONCAVE SPHERES line, it is advisable to replace the RETENTIVE CAP every 12 months.

Any use of the TITANIUM+TiN CONCAVE SPHERES line which does not follow the present instructions or other RHEIN 83 S.R.L. literature is considered improper.

Legend of the symbols:

	"Consult instructions for use"
	"Do not re-use"
	"Batch code"
	"Catalog number"
	"Manufacturer"
	"Non-sterile"
	"Keep away from sunlight"
	"Do not use if package is damaged"