Fabbricante: RHEIN 83 S.R.L. Via E. Zago, 10 A,B,C
40128 Bologna, Italia • Tel. 051-244510 • Fax 051-245238 •
N. Verde 800-901172

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Instruction for use PIVOT OT EQUATOR, PIVOT FLEX, PIVOT BLOCK, OT REVERSE 3 and CASTABLE PIVOTS

DEVICE DESCRIPTION AND INDICATIONS FOR USE: the PIVOT and OT REVERSE 3 line is composed by titanium devices; the PIVOT OT EQUATOR and FLEX have the retentive head with TiN coating as well as the OT REVERSE 3 pivot. The line of CASTABLE PIVOTS is composed by plastic devices for the root canal impression or for the construction of the Richmond crowns. **IMPORTANT:** these instructions for use are the most recent available; please read and store them carefully. **CONTRAINDICATIONS:** at the present state of knowledge there are no contraindications to the use of the present line of products; at the date of emission of these instructions for use no undesirable side effects have been reported. **CAUTION:** the choice of the attachment is demanded to the dentist or to the dental technician according to the prosthetic project. **MATERIAL:** all the materials used for the present line of products are of quality and meet the international standards of reference. The patient has to inform the dentist about potential allergies. If allergic reactions occur search the source and remove it. **WARNINGS AND SINGLE USE:** the present line of products is composed by single-use products; their reuse will undermine the correct functionality and may affect the patient's health. RHEIN 83 S.R.L. assumes no responsibility for any damage derived from an improper use or reuse of the present line of products. **TRANSPORTATION, STORAGE AND EXPIRATION DATE:** do not damage the packaging during the transportation; store in a clean, dry place, in the original packaging away from sunlight and heat sources. The product does not expire. **WASHING AND STERILIZATION:** the content of the set are sold as NON STERILE; the metal products should be washed in hot water with a soft brush and non-aggressive and non-corrosive detergent solutions (avoid solutions containing chlorine). The dentist is required to pack and store the prefabricated metal parts in cleared FDA 510(k) pouches (paper/plastic film) suitable for sterilizing and have to sterilize them before use with the following parameters: steam sterilization with prevacuum cycles at 135 °C – 3 minutes, drying time 16 minutes, pressure 2.06 bar. It is recommended to sanitize the retentive CAPS with benzalkonium chloride solutions for 10 minutes, to rinse with distilled water and dry them before use. The dentist is responsible of the sterilization. The castable components do not require sterilization as they are subjected to transformation (melting) by the dental technician. **TECHNICAL SUPPORT:** our technical staff is available for technical questions to the numbers listed above or to the email: gianni.storni@rhein83.it. For more details on the use of the present line of products, please consult the catalog-manual and the leaflets. **MAINTENANCE AND PERIODIC CARE:** dentists have the responsibility to keep the proper functionality and retention of the attachments assuring the safety of the patient by constant maintenance. In order to maintain the high quality standard offered by the present line of products we recommend the substitution of the retentive elastic components yearly. **GUIDELINES FOR THE PATIENTS:** patients are recommended to follow the indications provided by the dentist, to attend periodical controls and perform daily accurate hygiene. The present line of products have not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the present line of products in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. **DISPOSAL:** follow the local laws for the disposal. It is suggested to wash and sterilize the products accordingly to the surgical protocols before disposal. **CONTENT OF THE PACKAGES:** all the components of the present product line can be sold individually, only some of them are also part of sets.

TECHNICAL SPECIFICATIONS

PIVOT OT EQUATOR titanium pivot with 2.5 mm hemi-sphere available in various lengths of the conical body 7 mm, 9 mm (fuchsia), 10 mm (green): SINGLE USE titanium retentive posts for definitive use. The conical part is sandblasted and with retentions for definitive cementation with definitive cement; the pivots are anodized in green and fuchsia colors for a better identification of the various lengths. The various lengths of the conical part have the same conicity of the common burs for the preparation of the root canal. The retentive semi-sphere has a nitride coating and has an internal thread for the screwing of the OT EQUATOR IMPRESSION COPING with screw.

PIVOT BLOCK with Normal size sphere of 2.5mm and Micro size sphere of 1.8mm available in various lengths of conical body 7 mm, 9 mm (fuchsia), 10 mm (green): SINGLE USE titanium retentive posts for temporary or permanent use. The conical part is smooth for temporary cementation with provisional cement; the pivots are anodized in green and fuchsia colors for a better identification of the various lengths. The various lengths of the conical part have the same conicity of the common burs for the preparation of the root canal. For final cementation, follow the instructions in the paragraph "Permanent cementation of PIVOT BLOCK".

PIVOT FLEX mobile Normal size sphere of 2.5 mm available in various lengths of the conical body 7 mm, 9 mm (fuchsia), 10 mm (green): SINGLE USE titanium retentive posts for definitive use. The conical part is sandblasted and with retentions for definitive cementation with definitive cement. The various lengths of the conical part have the same conicity of the common burs for the preparation of the root canal. The retentive sphere has a nitride coating and can tilt up to 7.5° to overcome a greater disparallelism of the root.

OT REVERSE 3: SINGLE USE titanium TiN coated permanent retentive posts. The Micro RETENTIVE MALE in elastic material is provided with a titanium insert that provides retention to the denture. The STAINLESS STEEL HOUSING, will host the RETENTIVE MALE.

OT REVERSE 3 ACCESSORIES: IMPRESSION COPING, MANUAL TOOL, OT REVERSE 3 BUR and ANALOG for laboratory model.

PROTECTIVE DISK: SINGLE USE plastic devices available in two variants: transparent and pink; used to protect the gum and eliminate the undercuts during fixation with self-curing resin, of the metallic housing + retentive CAP in the mouth of the patient.

CASTABLE PIVOTS available in various lengths and diameter of the conical body: SINGLE USE castable plastic posts identified by different colors and different diameters of the spheres. Their different length allows a more accurate impression of the root canal.

INSTRUCTION FOR USE

The roots should be first treated with endodontical treatment and prepared to receive the PIVOT to be cemented or the castable PIVOT for the impression.

Temporary cementation of PIVOT BLOCK (Normal and Micro size): once the canal site is ready, choose the PIVOT with the correct length; it must be passive, but not loose, into the root canal. Insert the GRAY DIRECTIONAL RING on the PIVOT sphere, taking care to first insert the widest hole in order to properly seat the RING. Insert the IMPRESSION COPING with the CAP which is coupled to the GRAY RING allows you to easily handle the post in the cementation phase. In addition, the RING and the IMPRESSION COPING will avoid the cement to touch the surface of the sphere. Insert the temporary cement into the root canal and insert the post until the RING touches the root. Wait for the complete hardening of the temporary cement by carefully following the manufacturer's instructions. Once the temporary cement is hardened, remove the IMPRESSION COPING and the RING by removing any excess cement.

Permanent cementation of PIVOT BLOCK (Normal and Micro size) and PIVOT FLEX: once the canal site is ready, choose the PIVOT with the correct length; it must be passive, but not loose, into the root canal. Insert the GRAY DIRECTIONAL RING on the PIVOT sphere, taking care to first insert the widest hole in order to properly seat the RING. Then insert the IMPRESSION COPING with the CAP which is coupled to the GRAY RING allows you to easily handle the post in the cementation phase. In addition, the RING and IMPRESSION COPING will avoid the cement to touch the surface of the sphere. Insert the definitive cement into the root canal and insert the post until the RING touches the root. Wait for the complete hardening of the definitive cement by carefully following the manufacturer's instructions. Once the definitive cement is hardened, remove the IMPRESSION COPING and the RING by removing any excess cement.

Permanent cementation of PIVOT OT EQUATOR: once the canal site is ready, choose the PIVOT with the correct length; it must be passive, but not loose, into the root canal. Insert the IMPRESSION COPING with the CAP on the hemi-sphere that allows you to easily handle the post in the cementation phase. In addition, the IMPRESSION COPING will avoid the cement to touch the surface of the sphere. Insert the definitive cement into the root canal and insert the post until the IMPRESSION COPING touches the root. Wait for the complete hardening of the definitive cement by carefully following the manufacturer's instructions. Once the definitive cement is hardened, remove the IMPRESSION COPING and any excess of cement.

Fixing of the CAPS into the mouth (for PIVOT BLOCK and FLEX): Choose the RING with the inclination that is best suited to compensate the disparallelism of the attachment. Place it on the neck of the sphere. Pay attention to first insert the wider hole of the DIRECTIONAL RING; the narrowest hole will then face the sphere. Insert the PROTECTIVE DISC on the RING and then the HOUSING + CAP; turn the RING until the CAP will be perpendicular to the insertion plane or parallel to the occlusal plan, check that there is enough space for the HOUSING in the denture and cure the resin. If you proceed with the impression, insert the IMPRESSION COPING directly above the DIRECTIONAL RING and proceed with the impression.

Fixing of the CAPS into the mouth (for PIVOT OT EQUATOR): insert the PROTECTIVE DISC and then the HOUSING + RETENTIVE CAP; alternatively insert the SMART BOX housing with the SMART BOX black cap, check that there is enough space for the HOUSING in the denture and cure the resin. If you proceed with the impression, insert the IMPRESSION COPING directly above the PIVOT and proceed with the impression.

Permanent cementation of OT REVERSE 3 PIVOT and fixing of the RETENTIVE MALE: once the canal site is ready, ream the root canal with the OT REVERSE 3 BUR and check that the post is passive, it must be passive, but not loose, into the root canal. Insert the MANUAL TOOL that allows to comfortably manage the post during the cementing phase, avoiding that the cement will enter the hole of the post. Insert the permanent cement into the root canal and insert the post until the MANUAL TOOL touches the root acting as a stop. Wait for the complete hardening of the permanent cement strictly following the manufacturer's instructions. Once the permanent cement has hardened, remove the MANUAL TOOL and remove any excess cement. Insert the RETENTIVE MALE in the INOX Micro HOUSING, check that there is enough space in the denture for the HOUSING and cure the resin. If you proceed with taking the impression, insert the IMPRESSION COPING directly over the PIVOT and proceed with taking the impression.

CASTABLE PIVOTS: once the canal site is ready, choose the PIVOT with the correct length; it must be passive, but not loose, into the root canal. Brush the compatible adhesive to the selected impression material, put the impression material on the castable PIVOT and proceed with the impression. Once the working model is developed, use the right length of the CASTABLE PIVOT to create the post of the Richmond crown. Insulate the model channel and reline the pivot with the castable resin by inserting it into the canal. Once the resin is hardened complete the modeling of the Richmond crown by placing perpendicular to the occlusal plane the chosen castable attachment, sprue and cast. The dental technician is demanded for the choice of the metal (we recommend an alloy with hardness of at least 280 Vickers) and the investment material to be used in accordance with the operating instructions supplied by the respective manufacturers.

To avoid the loss of the performances of the present line of products, replace the retentive CAPS every 12 months.

Any use of the present line of products which does not follow the present instructions or other RHEIN 83 S.R.L. literature is considered improper.

Legend of the symbols:

	"Consult instructions for use"
	"Do not re-use"
	"Batch code"
	"Catalog number"
	"Manufacturer"
	"Non-sterile"
	"Keep away from sunlight"
	"Do not use if package is damaged"