RHEIN 83 S.R.L.

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Type: TITANIUM DIRECT ATTACHMENTS

Modello/Model: SPHERO FLEX – SPHERO BLOCK NORMO – SPHERO BLOCK
MICRO – OT EQUATOR – OT EQUATOR & SMARTBOX

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Rhein 83 S.R.L., in accordance with the new guidelines on SSCP management in EUDAMED as set forth in MDCG 2026-4 (Position Paper on the Management of SS(C)Ps in EUDAMED), is revising this document in the areas outlined in the guideline. The revisions are highlighted in green for ease of reading.

The changes/additions introduced by the guideline concern:

- Responsibility for uploading the SSCP—in the “master” version validated by the Notified Body and in the relevant translations pursuant to Article 29(4) of the MDR—which will be transferred from the Notified Body to the manufacturer.
- Translations and readability: this document is made available in the languages required by the Member States in which the devices are marketed, in accordance with the provisions for the Instructions for Use (IFU), and is in any case available in English

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Separate documents will be created for USERS and PATIENTS, which are compliant with and fully consistent with this document in terms of content, and which will be included on the respective company web pages.

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This Summary of Safety and Clinical Performance (SSCP) document is intended to provide public access to an up-to-date summary of the key safety and clinical performance aspects of the Medical Devices (herein M.D.) listed below.

The information presented below is intended for patients.

A more extensive summary prepared for healthcare professionals can be found in the first part of this document.

The SSCP document is not intended to provide general advice on the treatment of a medical condition.

Please contact your physician if you have any questions about your medical condition, the use of the device in your situation or the evaluation of possible alternative treatments.

This document is not intended to replace the Implant Passport (Implant Card) or the Instructions for Use with information on the safe use of the device.

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FIRST YEAR OF ISSUE OF THE EC CERTIFICATE MDR:
2023

NOTIFIED BODY WHO VALIDATE THE SSCP DOCUMENT AND ITS IDENTIFICATION NUMBER:
ICIM 0425

MANUFACTURER IDENTIFICATION:

RHEIN 83 S.R.L.

Via Zago 10/ABC

40128 Bologna (BO) – Italy

Tel.: 051 244510 – <https://www.rhein83.com> - info@rhein83.com

Single Registration Number (SRN): IT-MF-000005998

RHEIN 83 S.R.L. is responsible for placing on the EU market the M.D.s subject of this document.

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The M.D.s covered by this document consist of: titanium direct attachments with nitride coating. Designed to allow the practitioner to create a removable mucous-supported, implant-retained denture.

The Custom-made Medical Device (C.M.D.) replaces a missing anatomical structure to restore masticatory function, correct phonation, support of the perioral tissues, and last but not least, the smile and consequent relational ability.

NOMENCLATURA EUDAMED MDR - CND e UDI-DI di base/ NOMENCLATURE EUDAMED MDR – CND and Basic UDI-DI

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DM / MD	CLASSE / CLASS	REGOLA / RULE	UDI-DI di base / Basic UDI-DI
ATTACCO SPHERO BLOCK NORMO SPHERO BLOCK NORMO ATTACHMENT	IIB	8	805804571ADTK5
ATTACCO SPHERO BLOCK MICRO SPHERO BLOCK MICRO ATTACHMENT			
ATTACCO SPHERO FLEX SPHERO FLEX ATTACHMENT			
ATTACCO OT EQUATOR OT EQUATOR ATTACHMENT			
ATTACCO OT EQUATOR & SMARTBOX OT EQUATOR ATTACHMENT & SMARTBOX			

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These MDs are indicated for partially or fully edentulous patients who wear or do not wear dentures when the creation of a removable, implant-supported denture is prescribed.

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The manufacturer places no restrictions in terms of the selection of patients on whom DM products can be used as all products are biocompatible.

They can be used on patients of any age, predominantly adults or elderly (> 16-18 years) not suffering from uncontrolled systemic diseases or general surgical contraindications. Edentulous or partially edentulous adult patients are the target patients.

In rare cases, as a result of dental malformations or agenesis, it is possible for these patients to be developmentally disabled, as documented in Art. 103 "oral rehabilitation with implants..." of the scientific bibliography RHEIN 83 S.R.L. in which a case of prosthetic rehabilitation of a child aged 11 years with hypohidrotic ectodermal dysplasia (HED), a genetic defect of ectodermal development with malformations of the ectodermal structures (skin, hair, teeth and sweat glands) was documented. This case was developed at the Odontostomatological Sciences Unit of the University of Bologna.

The patient selection criteria are defined by the dentist/dentist team, which assesses the best prosthetic treatment plan based on the patient's medical history.

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These MDs are contraindicated for partially or totally edentulous patients whom the dentist/dentist team does not identify as suitable for the following prosthetic treatments: Mucus-supported removable denture and retained implant.

These DMs must be used by qualified personnel with specific dental or dental professional training. The prosthetic design must meet the specific function, aesthetic and phonetic needs of the intended patient.

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The operating principle and mode of action of these DMs is simple, known to the state of the art and established in use since 1983, when they were first marketed.


Regardless of the geometric shape of the sphere (SPHERO BLOCK NORMO & MICRO and SPHERO FLEX) or the hemisphere (OT EQUATOR), the operation is identical for all. In fact, these DMs act by exploiting the retention offered by the elastic cap that 'anchors' in the undercut created by the head of the direct attachment screwed into the dental implant. In other words, the retention is given by the effective contact surface between the apical part of the sphere, or semi-sphere, (undercut) and the elastic cap retained inside its metal container, which in turn is incorporated into the removable denture. In addition, to obtain different degrees of retention, copings of different elasticity identified by colour are produced.

Each cap is made by a specific material that gives it a greater or lesser retention. A careful study of materials was therefore carried out in order to choose products that could be biocompatible and offer safety and performance results appropriate to their intended use.


Elastic materials allow the cap to allow the ball or hemisphere to enter, and then return to its actual shape, thus achieving the wrap-around effect that allows for maximum retention.

The more than 30 years of use in dentistry and the diffusion in over 100 countries worldwide of these DMs guarantee excellent evidence of their effectiveness and safety. In fact, these DMs were created in response to the needs of dentist and thus from requirements 'from the field'. Their validity is supported by years of use by professionals all over the world who always have the opportunity to report any product nonconformities to the company. This information is collected by the company, together with feedback received during after-sales service (technical assistance from the R&D lab) to which feedback received from end-customer training courses (training courses held by the R&D lab and Opinion Leaders) is added. All this information is periodically evaluated through the various post-market surveillance and risk reassessment actions.

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

OT EQUATOR scientifically proven mode of action:

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SPHERO BLOCK NORMO, MICRO e FLEX scientifically proven mode of action:

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The devices in question are legacy devices, already certified according to Directive 93/42/EEC by the notified body ICIM with certificate No. 0425-MED-003705-00 and already on the EU and international markets.

No significant changes have been made to the devices compared to the devices already certified according to the Directive.

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Disposable plastic accessory. **The directional rings (grey ring 0°, green ring 7°, red ring 14°) are invasive of the oral cavity for the time strictly necessary when the housing with cap is fixed inside the removable denture, in the patient's mouth.**

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
SINGLE-USE plastic accessory. Black cap-positioner designed to align and protect from liquid resin infiltration during the direct embodying phase of the Smartbox housing inside the denture. It is therefore invasive of the oral cavity for the time strictly necessary for this operation which, once completed, will see the replacement of the black cap-positioner with an OT Equator retentive cap.

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Nitrided titanium device STERILIZABLE in an autoclave. **Hexagonal connection driver for screwing/unscrewing Sphero Block (Normo & Micro) and Sphero Flex attachments** in the mouth or on the laboratory model. Compatible with RHEIN 83 dynamometric ratchet or dynamometric handpiece connector.

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Reversible manual torque wrench. **For tightening retentive abutments or implant prosthetic screws to a predetermined torque.**

Can also be customized with the laser marked customer's logo.

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SQUARE KEY WITH HOLDER FOR OT EQUATOR CL.I R.5

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REUSABLE DM made of stainless steel. The screwdriver has been heat-treated and is autoclavable; available in long and short versions with a square-conical tip and a rotating plastic ring that facilitates the tightening process. Compatible with the RHEIN 83 S.R.L. torque ratchet. **driver is used to screw/unscrew the OT Equator attachments** in the mouth or on the laboratory model, manually or by means of a torque wrench.

CONNETTORE QUADRATO-CONICO PER MANIPOLO DINAMOMETRICO CL.IIA R.6
SQUARE-CONIC SCREWDRIVER FOR DYNAMOMETRIC HANDPIECE CL.IIA R.6

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Reusable DM made of stainless steel. The screwdriver has been heat-treated and is autoclavable; available in long and short versions with a square-conical tip. **The SCREWDRIVER is used to screw/unscrew the OT Equator attachments** in the mouth or on the laboratory model, using a dynamometric handpiece.

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REUSABLE and STERILIZABLE titanium device with lacquered reference marks **to measure the height of the cuff height in order to choose the correct abutment height.**

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They have not been evaluated for safety and compatibility in the MRI environment.

They have not been tested for heating, migration, or imaging artefacts in the MRI environment. The safety of these M.D.s in the MRI environment is unknown.

Scanning a patient with these M.D.s could result in patient injury.

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The dentist must check for allergies before starting prosthetic rehabilitation. In the case of allergic events, search for the cause and remove it. In the event of a proven allergy to the materials constituting these M.D.s, do not use these M.D.s. The dentist should inform the patient of any new M.D.-related event; if a serious incident occurs, this should be reported immediately to the Competent Authority and the Manufacturer.

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The dentist must provide the intended patients with instructions for correct daily hygiene and appropriate use of the D.M. and D.M.M. Do not use these devices on patients with an intolerance or allergy to titanium. Users must carry out the prosthesis wearing all the protective equipment provided for their profession (e.g., mask, gowns, gloves, goggles) and take the utmost care to avoid minor injuries from contact with devices with sharp or pointed portions and to avoid contaminating themselves with the patient's saliva or blood.

RHEIN 83 S.R.L. provides dental practitioners with an Implant Passport (Implant Card) which, after being properly filled out and completed with the labels of the implants inserted, must be given to the patient.

In this document all information regarding the implant is recorded, which is useful for the traceability of the device over time.

If necessary, the patient can show the Implant Passport to any dentist who will then be able to intervene appropriately in resolving any problems and contact the manufacturer if necessary.

It is suggested that the dentist also make a note on the Implant Passport of the type of prosthetic rehabilitation performed, the date on which it was performed and if possible, the codes and lots of the prosthetic components used.

It is the patient's duty to carefully keep and carry the Implant **Passport with him/her at each subsequent check-up.**