RHEIN 83 S.R.L.

FASCICOLO TECNICO 02 MDR / MDR TECHNICAL FILE 02

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Type: TITANIUM DIRECT ATTACHMENTS

Modello/Model: SPHERO FLEX – SPHERO BLOCK NORMO – SPHERO BLOCK
MICRO – OT EQUATOR – OT EQUATOR & SMARTBOX

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Rhein 83 S.R.L., in accordance with the new guidelines on SSCP management in EUDAMED as set forth in MDCG 2026-4 (Position Paper on the Management of SS(C)Ps in EUDAMED), is revising this document in the areas outlined in the guideline. The revisions are highlighted in green for ease of reading.

The changes/additions introduced by the guideline concern:

- Responsibility for uploading the SSCP—in the “master” version validated by the Notified Body and in the relevant translations pursuant to Article 29(4) of the MDR—which will be transferred from the Notified Body to the manufacturer.
- Translations and readability: this document is made available in the languages required by the Member States in which the devices are marketed, in accordance with the provisions for the Instructions for Use (IFU), and is in any case available in English

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Separate documents will be created for USERS and PATIENTS, which are compliant with and fully consistent with this document in terms of content, and which will be included on the respective company web pages.

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This Summary of Safety and Clinical Performance (SSCP) document is intended to provide public access to an up-to-date summary of the key safety and clinical performance aspects of the Medical Devices (herein M.D.) listed below. The information presented below is intended for professional users.

FIRST YEAR OF ISSUE OF THE EC CERTIFICATE MDR:
2023

NOTIFIED BODY WHO VALIDATE THE SSCP DOCUMENT AND ITS IDENTIFICATION NUMBER:
ICIM 0425

MANUFACTURER IDENTIFICATION:

 **RHEIN 83 S.R.L.**
Via Zago 10/ABC

40128 Bologna (BO) – Italy

Tel.: 051 244510 – <https://www.rhein83.com> - info@rhein83.com

Single Registration Number (SRN): IT-MF-000005998

[RHEIN 83 S.R.L.](#) is responsible for placing on the EU market the M.D.s subject of this document.

MEDICAL DEVICE IDENTIFICATION:

SPHERO FLEX - SPHERO BLOCK NORMO - SPHERO BLOCK MICRO - OT EQUATOR - OT EQUATOR & SMARTBOX

The M. D.s covered by this document consist of: titanium direct attachments with nitride coating. Designed to allow the practitioner to create a removable mucous-supported, implant-retained denture.

The Custom-made Medical Device (C.M.D.) replaces a missing anatomical structure to restore masticatory function, correct phonation, support of the perioral tissues, and last but not least, the smile and consequent relational ability.

NOMENCLATURE EUDAMED MDR – CND and Basic UDI-DI

MDN 1103 (NON-ACTIVE DENTAL IMPLANTS AND DENTAL MATERIAL)

P010299 (ODONTOLOGICAL PROTHESES – OTHER)

NOMENCLATURE EUDAMED MDR – CND and Basic UDI-DI

MDN 1103 (NON-ACTIVE DENTAL IMPLANTS AND DENTAL MATERIAL)

P010299 (ODONTOLOGICAL PROTHESES – OTHER)

DM / MD	CLASSE / CLASS	REGOLA / RULE	UDI-DI di base / Basic UDI-DI
ATTACCO SPHERO BLOCK NORMO SPHERO BLOCK NORMO ATTACHMENT	IIB	8	805804571ADTK5
ATTACCO SPHERO BLOCK MICRO SPHERO BLOCK MICRO ATTACHMENT			
ATTACCO SPHERO FLEX SPHERO FLEX ATTACHMENT			
ATTACCO OT EQUATOR OT EQUATOR ATTACHMENT			
ATTACCO OT EQUATOR & SMARTBOX OT EQUATOR ATTACHMENT & SMARTBOX			

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INTENDED USE

These M.D.s have been designed to allow the practitioner to create a removable, mucus-supported, implant-retained denture.

The C.M.D. replaces a missing anatomical structure to restore masticatory function, correct phonation, support of the perioral tissues, and last but not least, the smile and consequent ability to relate.

This definition is consistent with what is recited in point 1 of Article 2 "Definitions" in CHAPTER I "Scope and Definitions" of the MDR EU 2017/745, which defines (in bold the parts applicable to these MDRs):

"medical device" means any instrument, apparatus, appliance, software, implant, reagent, material or other article, intended by the manufacturer to be used on humans, alone or in combination, for one or more of the following specific medical purposes:

- diagnosis, prevention, monitoring, prediction, prognosis, treatment or mitigation of disease,
- diagnosis, monitoring, treatment, alleviation or compensation for an injury or disability,
- study, replacement or modification of anatomy or of a physiological or pathological process or state,
- providing information through the in vitro examination of samples from the human body, including donated blood and tissue,

and which does not exert its principal intended action in or on the human body by pharmacological, immunological or metabolic means, but whose function may be assisted by such means.

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CLASSIFICATION

M.D. Type	Description	Single Use	Class	Images
<p>SPHERO BLOCK NORMO ATTACHMENT</p>	<p>Single-use DM consisting of:</p> <ul style="list-style-type: none"> -a nitride titanium direct attachment with a Ø2.5mm retention ball; -a series of plastic caps, with an internal diameter of 2.5 mm, which are anchored to the retention ball. The attachment is screwed into a dental implant using a special key; - stainless steel housing; -a series of directional rings (grey, green and red) for the correct positioning of the metal container + cap inside the denture; -a protective disc to avoid contact between the liquid resin and the gingiva during the phase of fixing the metal container + coping inside the denture. <p>The denture can be anchored onto the implants with disparallelism of to 14° per implant. This guarantees adequate retention in the patient's mouth.</p>	<p>YES</p>	<p>II B</p>	
<p>SPHERO BLOCK MICRO ATTACHMENT</p>	<p>Single-use DM consisting of:</p> <ul style="list-style-type: none"> -a nitride titanium direct attachment with a Ø1.8mm retention ball; -a series of plastic caps, with an internal diameter of 1.8 mm, which are anchored to the retention ball. The attachment is screwed into a dental implant using a special key; - stainless steel housing; -a series of directional rings (grey, green and red) for the correct positioning of the metal container + cap inside the denture; -a protective disc to avoid contact between the 	<p>YES</p>	<p>II B</p>	

	<p>liquid resin and the gingiva during the phase of fixing the metal container + coping inside the denture.</p> <p>The denture can be anchored onto the implants with disparallelism of to 14° per implant. This guarantees adequate retention in the patient's mouth.</p>			
<p>SPHERO FLEX ATTACHMENT</p>	<p>Single-use DM consisting of:</p> <ul style="list-style-type: none"> - a direct nitride titanium attachment with a \varnothing 2.5 mm retention ball; this ball has the ability to rotate around the main axis of the attachment for 7.5° in every direction, being guaranteed the impossibility of detachment from the rotation pin; this mechanism facilitates the insertion of the denture on the attachments, even in the case of implants positioned with divergent inclinations. The attachment is screwed into a dental implant using a special key; - a series of plastic caps, with an internal diameter of 2.5 mm, which are anchored to the retention ball; - stainless steel housing; - a set of directional rings (grey, green and red) for the correct positioning of the metal container + cap inside the denture; - a protective disc to avoid contact between the liquid resin and the gingiva during the phase of fixing the metal container + coping inside the denture. <p>The denture can be anchored onto the implants with disparallelism of to 21,5° per implant. This guarantees adequate retention in the patient's mouth.</p>	<p>YES</p>	<p>II B</p>	
<p>OT EQUATOR ATTACHMENT</p>	<p>Single-use DM consisting of:</p> <ul style="list-style-type: none"> - a direct nitride titanium attachment with a \varnothing 2.5 mm retention hemisphere; the attachment is screwed into a dental implant using a special key; - a series of plastic caps, with an internal diameter of 2.5 mm, which are anchored to the retention hemisphere; - stainless steel housing; - a protective disc to prevent contact between the liquid resin and the gingiva during the phase of fixing the metal container + coping inside the denture. <p>The denture can be anchored onto the implants with disparallelism of to 15° per implant. This</p>	<p>YES</p>	<p>II B</p>	

	guarantees adequate retention in the patient's mouth.			
OT EQUATOR ABUTMENT & SMARTBOX	<p>Single-use DM consisting of:</p> <ul style="list-style-type: none"> - a direct nitride titanium attachment with a \varnothing 2.5 mm retention hemisphere; the attachment is screwed into a dental implant using a special key; - a series of plastic caps, with an internal diameter of 2.5 mm, which are anchored to the retention hemisphere; - titanium tilting housing; - a protective disc to prevent contact between the liquid resin and the gingiva during the phase of fixing the metal container + coping inside the denture. <p>The denture can be anchored onto the implants with disparallelism of to 25° per implant. This guarantees adequate retention in the patient's mouth.</p>	YES	II B	

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These M.D.s are indicated for partially or fully edentulous patients who wear or do not wear dentures when the creation of a removable, implant-supported denture is prescribed.

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The manufacturer places no restrictions in terms of the selection of patients on whom DM products can be used as all products are biocompatible.

They can be used on patients of any age, predominantly adults or elderly (> 16-18 years) not suffering from uncontrolled systemic diseases or general surgical contraindications. Edentulous or partially edentulous adult patients are the target patients.

In rare cases, as a result of dental malformations or agenesis, it is possible for these patients to be developmentally disabled, as documented in Art. 103 "oral rehabilitation with implants..." of the scientific

bibliography RHEIN 83 S.R.L. in which a case of prosthetic rehabilitation of a child aged 11 years with hypohidrotic ectodermal dysplasia (HED), a genetic defect of ectodermal development with malformations of the ectodermal structures (skin, hair, teeth and sweat glands) was documented. This case was developed at the Odontostomatological Sciences Unit of the University of Bologna.

The patient selection criteria are defined by the dentist/dentist team, which assesses the best prosthetic treatment plan based on the patient's medical history.

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These M.D.s are contraindicated for partially or totally edentulous patients whom the dentist/dentist team does not identify as suitable for the following prosthetic treatments: Mucus-supported removable denture and retained implant.

These DMs must be used by qualified personnel with specific dental or dental professional training. The prosthetic design must meet the specific function, aesthetic and phonetic needs of the intended patient.

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The operating principle and mode of action of these DMs is simple, known to the state of the art and established in use since 1983, when they were first marketed.

Regardless of the geometric shape of the sphere (SPHERO BLOCK NORMO & MICRO and SPHERO FLEX) or the hemisphere (OT EQUATOR), the operation is identical for all. In fact, these DMs act by exploiting the retention offered by the elastic cap that 'anchors' in the undercut created by the head of the direct attachment screwed into the dental implant. In other words, the retention is given by the effective contact surface between the apical part of the sphere, or semi-sphere, (undercut) and the elastic cap retained inside its metal container, which in turn is incorporated into the removable denture. In addition, to obtain different degrees of retention, copings of different elasticity identified by colour are produced.

Each cap is made by a specific material that gives it a greater or lesser retention. A careful study of materials was therefore carried out in order to choose products that could be biocompatible and offer safety and performance results appropriate to their intended use.

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



Elastic materials allow the cap to allow the ball or hemisphere to enter, and then return to its actual shape, thus achieving the wrap-around effect that allows for maximum retention.

The more than 30 years of use in dentistry and the diffusion in over 100 countries worldwide of these DMs guarantee excellent evidence of their effectiveness and safety. In fact, these DMs were created in response to the needs of dentist and thus from requirements 'from the field'. Their validity is supported by years of use by professionals all over the world who always have the opportunity to report any product nonconformities to the company. This information is collected by the company, together with feedback received during after-sales service (technical assistance from the R&D lab) to which feedback received from end-customer training courses (training courses held by the R&D lab and Opinion Leaders) is added. All this information is periodically evaluated through the various post-market surveillance and risk reassessment actions.

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

OT EQUATOR scientifically proven mode of action:

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SPHERO BLOCK NORMO, MICRO e FLEX scientifically proven mode of action:

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REFERENCE TO PREVIOUS GENERATIONS OR VARIANTS AND DESCRIPTION OF DIFFERENCES

The devices in question are legacy devices, already certified according to Directive 93/42/EEC by the notified body ICIM with certificate No. 0425-MED-003705-00 and already on the EU and international markets. No significant changes have been made to the devices compared to the devices already certified according to the Directive.

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Reversible manual torque wrench. For tightening retentive abutments or implant prosthetic screws to a predetermined torque.

Can also be customized with the laser marked customer's logo.

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SQUARE KEY WITH HOLDER FOR OT EQUATOR CL.I R.5

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REUSABLE DM made of stainless steel. The screwdriver has been heat-treated and is autoclavable; available in long and short versions with a square-conical tip and a rotating plastic ring that facilitates the tightening process. Compatible with the RHEIN 83 S.R.L. torque ratchet. **driver is used to screw/unscrew the OT Equator attachments** in the mouth or on the laboratory model, manually or by means of a torque wrench.

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SQUARE-CONIC SCREWDRIVER FOR DYNAMOMETRIC HANDPIECE CL.IIA R.6

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Reusable DM made of stainless steel. The screwdriver has been heat-treated and is autoclavable; available in long and short versions with a square-conical tip. **The SCREWDRIVER is used to screw/unscrew the OT Equator attachments** in the mouth or on the laboratory model, using a dynamometric handpiece.

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They have not been evaluated for safety and compatibility in the MRI environment.

They have not been tested for heating, migration, or imaging artefacts in the MRI environment. The safety of these M.D.s in the MRI environment is unknown.

Scanning a patient with these M.D.s could result in patient injury.

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The dentist must check for allergies before starting prosthetic rehabilitation. In the case of allergic events, search for the cause and remove it. In the event of a proven allergy to the materials constituting these M.D.s, do not use these M.D.s.

The dentist should inform the patient of any new M.D.-related event; if a serious incident occurs, this should be reported immediately to the Competent Authority and the Manufacturer.

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The devices are not suitable where a totally rigid connection is required; the use of the TRADITIONAL CONTAINER on an implant with more than 15° divergence is not recommended; the use of the SMARTBOX CONTAINER on an implant with more than 25° divergence is not recommended.

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The dentist must provide the intended patients with instructions for correct daily hygiene and appropriate use of the D.M. and D.M.M. Do not use these devices on patients with an intolerance or allergy to titanium. Users must carry out the prosthesis wearing all the protective equipment provided for their profession (e.g., mask, gowns, gloves, goggles) and take the utmost care to avoid minor injuries from contact with devices with sharp or pointed portions and to avoid contaminating themselves with the patient's saliva or blood.

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The use and handling of M.D.s covered by this document are reserved for medical and dental personnel with the necessary professional qualifications and training. Prosthetic reconstruction procedures must be carried out in a suitable environment with suitable asepsis. It is recommended to always cover the surfaces with sterile drapes, to cover the dental unit and the micromotor with suitable coverings, to isolate the operating field by covering the patient with suitable gowns, to wear sterile gloves, to open the instruments from their sterile wrappings just before their use

VALIDAZIONE CLINICA E FOLLOW-UP CLINICO POST-MARKETING/ CLINICAL VALIDATION AND POST-MARKETING CLINICAL FOLLOW-UP

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The clinical evaluation is updated annually, with the integration of data from the PMCF activity or in case of changes to DMs, in case of introduction of new devices, regulatory changes, adverse events detected on RHEIN 83 S.R.L. M.D. or equivalent competing Devices.

The devices in question are devices that have already been on the market in the EU for several years and with good clinical results (in accordance with Directive 93/42/EEC) and meet the General Safety and Performance Requirements required by Regulation 2017/745 with regard to the risk/benefit trade-off aspects.

The analysis of the PMCF activity, the analysis of the literature carried out that allowed to define the state of the art of safety and expected performance for these devices for prosthetic reconstruction after the application of dental implants and the considerations that emerged from the risk analysis and post-market surveillance allow to consider these devices suitable for their intended use.

COERENZA CON PSUR, PMCF E DATI DI VIGILANZA EUDAMED/ CONSISTENCY WITH PSUR, PMCF AND EUDAMED VIGILANCE DATA

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The safety and performance data reported in this document are periodically checked for consistency against the most recent PSUR and PMCF Evaluation Report and against the vigilance reports available in EUDAMED, in order to ensure alignment between the SSCP and the entire technical documentation of the device. This check is performed at least at each annual update of the PSUR and PMCF, in accordance with Art. 86 MDR and MDCG 2019-9 Rev.1 guidance.

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Until the end of 2025:

- No adverse events or non-serious incidents have occurred for M.D. manufactured by RHEIN 83 S.R.L. belonging to the Overdenture category (OT Equator, Sphero Block, Sphero Flex)
- Conditions requiring mandatory information to the competent authorities of any non-serious incidents and related corrective actions did not occur.
- There were no non-serious incidents to be reported to the MDSAP authorities
- As in previous years, there were no cases of Recall from the Market due to accidents or non-accidents.

The RHEIN 83 S.R.L. company has always characterized its relationship with customers and users of its devices with a high level of exchange of information and feedback on its medical devices.

Up to the end of 2025, no negative feedbacks or feedbacks were collected either on social media or directly in the field that could require a revision of the risk analysis.

There have been no recalls from the market due to accidents or near misses.

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In general:

- the quality of RHEIN 83 S.R.L. products continues to be the aspect most appreciated by customers, as well as the work of the sales offices and the punctuality of order deliveries

- No particularly significant suggestions were collected during this customer satisfaction analysis

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Analyzing sales data for all M.D. marketed by RHEIN 83 S.R.L.:

- There is no evidence of a significant number of complaints in relation to the sales volume of the devices produced

- Side effects attributable to the use of any medical device are to be considered absent

The following table shows the data from the clinical follow-up analysis, divided by year, where the following are shown:

- Sales volume EU countries;
- Serious incident complaints to be reported to the relevant bodies;
- Serious side effects to be reported to the relevant bodies.

The safety performance is underlined by the number of serious incidents to be reported to the relevant bodies, in relation to the total number of devices sold to the sample of users over the years as shown in the table below.

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It is therefore considered that the design choices made are in line with the state of the art in the field of prosthetic reconstruction.

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The alternative to the removable "overdenture" prosthesis is the conventional total denture which, when it does not achieve a sufficient degree of retention in the mouth, is covered with adhesive paste by the patient who accepts, however, the inconvenience of less retention as well as the need to reapply adhesive pastes, which always have a limited duration and less effective retention.

In view of the above, with only two implants it is reasonable, in terms of cost/benefit, to opt for a removable overdenture type prosthesis with direct attachment.

An alternative of greater caliber is the fixed prosthesis, which can only be made when the number and position of the implants make it a safe solution in bearing masticatory loads.

NORME DI RIFERIMENTO, LINEE GUIDA E SPECIFICHE COMUNI/ REFERENCE STANDARDS, GUIDELINES AND COMMON SPECIFICATIONS

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The harmonized standards with respect to EU Regulation 2017/745, applied to the M.D. in question are:

NAME	TITLE
ISO 13485:2016 / A11:2021 (ITA + ENG)	Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 15223-1:2025	Medical devices - Symbols to be used in information to be provided by the manufacturer - Part 1: General requirements
EN ISO 14971:2022	Medical devices - Application of risk management to medical devices

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Other non-harmonised standards considered when demonstrating compliance with the requirements of Regulation 2017/745:

NAME	TITLE
ASTM F136 - 13	Standard specification for Wrought Titanium-6 Aluminium-4 Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401)
ASTM F899 - 20	Standard Specification for Wrought Stainless Steels for Surgical Instruments
UNI EN ISO 10993-1:2021	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
UNI EN ISO 14801:2017	Dentistry - Implants - Dynamic fatigue testing of endosseous dental implants
MDR 745/2017/UE	COMMISSION REGULATION (EU) 2021/2226 (14/12/2021) of 14 December 2021 implementing Regulation (EU) 2017/745 of the European Parliament and of the Council as regards electronic instructions for use of medical devices

Linee guida MDCG - Regolamento 2017/745/ MDCG Guidelines - Regulation 2017/745:

NAME	TITLE
MDCG 2018-1	Guidance on BASIC UDI-DI and changes to UDI-DI
MDCG 2018-7	Provisional considerations regarding language issues associated with the UDI
MDCG 2019-8 Rev.2	Implant Card relating to the application of Article 18 Regulation (EU) 2017/745
MDCG 2021-11	Guidance on Implant Card – ‘Device types’
MDCG 2021-19	Guidance note integration of the UDI within an organization’s quality
MDCG 2021-24	Guidance on classification of medical devices
MDCG 2026-4	Position Paper on the management of SS(C)Ps in EUDAMED – upload

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Common Specifications:

No applicable common specifications are available at the time of writing.