



T H E D E N T A L A T T A C H M E N T S C O M P A N Y



1983-2023

CODE LIST

FIXED AND REMOVABLE PROSTHESIS

for dentists and dental technicians



RHEIN83

FIXED PROSTHESIS OT BRIDGE

<i>OT EQUATOR abutment, Cuff height measurer, KIT 42, Titanium abutments</i>	Pag. 3
<i>Castable abutments, Titanium abutments with lateral holes, Titanium abutment for endoral welding</i>	Pag. 4
<i>Healing abutments, Impression copings</i>	Pag. 4
<i>Laboratory analogs, Scan Abutment, Screws, OT EQUATOR screwing tools</i>	Pag. 5
<i>Screwing tools, Universal handle, Seeger inserter, T-Bar Components and accessories</i>	Pag. 6
<i>Spare parts, Gingival connectors, Seeger seat cleaning cutter, Ot Cem</i>	Pag. 7

REMOVABLE PROSTHESIS

<i>OT EQUATOR for implants, OT EQUATOR for bars, OT EQUATOR caps kit, OT EQUATOR castable</i>	Pag. 8
<i>Spare parts OT EQUATOR</i>	Pag. 9
<i>OT EQUATOR transfer and analogs, OT EQUATOR tools and accessories</i>	Pag. 10
<i>Single threaded spheres, Spare parts for single threaded spheres</i>	Pag. 11
<i>Sphero block and Sphero flex pfor implants, Cuff height measurer</i>	Pag. 12
<i>Tools Sphero flex and Sphero block, Anti-unscrewing system</i>	Pag. 13
<i>Laboratory set</i>	Pag. 14
<i>Basic laboratory set</i>	Pag. 15
<i>OT Cap / OT Box attachments, OT Box large attachments, OT Strategy attachments, OT Bar attachments</i>	Pag. 16
<i>OT Vertical attachments, OT Unilateral attachments, OT Cap Tecno attachments, OT Lock attachments</i>	Pag. 17
<i>Pivots set for clinic, OT EQUATOR pivots set</i>	Pag. 18
<i>Pivot OT EQUATOR caps set, Pivot Block / Flex</i>	Pag. 19
<i>Spare parts OT Cap normal</i>	Pag. 20
<i>Spare parts OT Cap micro</i>	Pag. 21
<i>OT Cap assortment sets</i>	Pag. 22
<i>Tools and accessories OT Cap normal / micro</i>	Pag. 23
<i>OT Strategy spare parts</i>	Pag. 24-25
<i>OT Bar multiuse spare parts</i>	Pag. 25
<i>OT Vertical spare parts, OT Reverse spare parts</i>	Pag. 26
<i>OT Lock spare parts, Sphere repair sets</i>	Pag. 27
<i>OT Lever, Small parallelometer, Denture boxes</i>	Pag. 28
<i>OT Cem, Broken screw extractor kit</i>	Pag. 29
<i>Rhein83: The training</i>	Pag. 30
<i>Rhein83 direct services</i>	Pag. 31

OT EQUATOR ABUTMENT



Ref. 030

- 1 OT EQUATOR TITANIUM ABUTMENT WITH TiN COATING

CUFF HEIGHT MEASURER



Ref. 008MBG

SET CONTENT

- 1 CUFF HEIGHT SLIDER GAUGE
- 1 CUFF HEIGHT FIXED ROD GAUGE
- 1 SILICON RING DISPENSER
- 20 SILICON RINGS

KIT 42



Ref. BE42
+ n° 42 OT
EQUATOR
(Ref. 030 to be chosen)

- 1 ASSORTMENT SET (autoclavable) + 42 OT EQUATOR +TiN ABUTMENTS
THE SET42 CAN BE CUSTOMIZED ON THE BASIS OF THE IMPLANTS BRANDS, IMPLANTS
DIAMETERS AND CUFF HEIGHTS DESIRED. PLEASE INDICATE, IN THE ORDER, THE QUANTITY
AND THE ABOVE SPECIFICATIONS FOR EVERY KIND OF ABUTMENT.
(Tools not included in the set)



Ref. BE42

- 1 EMPTY BOX FOR ASSORTED OT EQUATOR ABUTMENTS

ABUTMENTS FOR FIXED PROSTHESIS



Ref. MTSEG

Ø3,5mm

Ref. MT40SEG

Ø4mm

Ref. MT50SEG

Ø5mm

EXTRAGRADE TITANIUM ABUTMENT

SET CONTENT

- 1 EXTRAGRADE TITANIUM ABUTMENT TO BE CEMENTED h9,5mm
- 1 TITANIUM SCREW h2mm
- 2 WHITE SEEGER WITH HANDLE
- 1 PINK SEEGER WITH HANDLE



Ref. MMSEG

Ø3,5mm

Ref. MM40SEG

Ø4mm

MONCONE IN TITANIO EXTRAGRADE MINI

SET CONTENT

- 1 MINI EXTRAGRADE TITANIUM ABUTMENT TO BE CEMENTED h5,7mm
- 1 TITANIUM SCREW FOR MINI EXTRAGRADE TITANIUM ABUTMENT
- 2 WHITE SEEGER WITH HANDLE
- 1 PINK SEEGER WITH HANDLE



Ref. MTSFEG

Ø3,5mm

Ref. MTEG15SF

Ø4mm

EXTRAGRADE TITANIUM ABUTMENT INCLINED AT 15°

SET CONTENT

- 1 EXTRAGRADE TITANIUM ABUTMENT INCLINED AT 15° TO BE CEMENTED
- 2 WHITE SEEGER WITH HANDLE
- 1 PINK SEEGER WITH HANDLE
- 1 SPHERICAL TITANIUM SCREW FOR 15° INCLINED ABUTMENT

ABUTMENTS FOR FIXED PROSTHESIS



Ref. CMEG40C

Ø4mm

EXTRAGRADE CASTABLE ABUTMENT

SET CONTENT

- 1 EXTRAGRADE CASTABLE ABUTMENT WITH HOLE
- 1 TITANIUM SCREW h2 mm
- 2 WHITE SEEGER WITH HANDLE
- 1 PINK SEEGER WITH HANDLE



Ref. MCEG15SF

Ø4mm

EXTRAGRADE CASTABLE ABUTMENT WITHOUT HOLE

SET CONTENT

- 1 EXTRAGRADE CASTABLE ABUTMENT WITHOUT HOLE (ONLY AVAILABLE WITH Ø4mm)
- 2 WHITE SEEGER WITH HANDLE
- 1 PINK SEEGER WITH HANDLE



Ref. CMTFP

h1,8mm

Ref. CMTB30FP

h3mm

TITANIUM ABUTMENT WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE WITH WIRE-REINFORCEMENT

SET CONTENT

- 1 TITANIUM ABUTMENT WITH LATERAL HOLES (ONLY AVAILABLE WITH Ø 3,5mm)
- 1 TITANIUM SCREW h2 mm
- 2 WHITE SEEGER WITH HANDLE
- 1 PINK SEEGER WITH HANDLE



Ref. MT40GL

Ø4mm

TITANIUM ABUTMENT FOR ENDORAL WELDING

SET CONTENT

- 1 TITANIUM NORMAL ABUTMENT FOR INTRA ORAL WELDING h9,5 mm Ø4 mm
- 1 TITANIUM SCREW h2 mm
- 2 WHITE SEEGER WITH HANDLE
- 1 PINK SEEGER WITH HANDLE

HEALING ABUTMENTS



Ref. 159VAG

SCREW AND HEALING RING

SET CONTENT

- 1 TITANIUM SCREW
- 1 TITANIUM HEALING RING



Ref. 159EAG

OT EQUATOR AND HEALING RING

SET CONTENT

- 1 TITANIUM + TiN OT EQUATOR SCREW FOR HEALING RING
- 1 TITANIUM HEALING RING



Ref. 159VEA

- 1 Titanium + TiN OT EQUATOR screw
(for self-extracting peek ring or titanium self-extracting healing ring)

IMPRESSION COPINGS



Ref. 144TLE

- 1 TITANIUM OT EQUATOR LONG IMPRESSION COPING + TITANIUM LONG SCREW h11,5mm



Ref. 144TTE

- 1 TITANIUM OT EQUATOR IMPRESSION COPING + TITANIUM SCREW h6 mm



Ref. 144MTE

- 2 OT EQUATOR IMPRESSION COPINGS FOR INDIVIDUAL TRAY



Ref. TPM2

- 2 MINI OT EQUATOR IMPRESSION COPINGS h3,7mm

LABORATORY ANALOGS



Ref. 144AE

- 2 LABORATORY ANALOGS ø 4mm

Ref. AI502

- 2 LABORATORY ANALOGS ø 5mm



Ref. 144AVC4

- 1 LABORATORY ANALOG WITH SCREW FOR CAD/CAM ø 4mm

Ref. 144AVC5

- 1 LABORATORY ANALOG WITH SCREW FOR CAD/CAM ø 5mm

OT EQUATOR SCAN ABUTMENT



Ref. 145SAE

- 1 OT EQUATOR TITANIUM SCAN ABUTMENT + TITANIUM SCREW

SCREWS



Ref. 147VME

- 1 TITANIUM SCREW FOR EXTRAGRADE ABUTMENT
(for screwdriver/stem for inclined Channel)



Ref. 146VMI

- 1 SPHERICAL TITANIUM SCREW FOR INCLINED ABUTMENT
(for screwdriver/stem for inclined Channel)



Ref. 148VCM

- 1 TITANIUM SCREW FOR MINI EXTRAGRADE ABUTMENT
(for screwdriver/stem for inclined Channel)



Ref. VC

- 1 TITANIUM SCREW FOR EXTRAGRADE ABUTMENT h2mm



Ref. VM

- 1 TITANIUM SCREW FOR EXTRAGRADE ABUTMENT h6mm



Ref. VXL

- 1 TITANIUM SCREW FOR EXTRAGRADE ABUTMENT h11,5mm

TOOLS



Ref. 774CHE

Ref. 775CLE*

(*Long version)

- 1 SQUARE SCREW DRIVER + HOLDER for tightening OT Equator
(SQUARE TIP 1.25MM) (USABLE WITH THE MANUAL TORQUE WRENCH)



Ref. 774HC

- 1 INTERCHANGEABLE HOLDER



Ref. 760CSC

Ref. 760CS*

(*Long version)

- 1 SCREWDRIVER/STEM FOR ANGULATED CHANNEL
FOR TORQUE CONTROLLER CONTRA ANGLE (SQUARE 1.25MM)
AND ADAPTER FOR MANUAL TORQUE WRENCH

TOOLS



Ref. 760ACC

- 1 SCREWDRIVER/STEM FOR ANGULATED CHANNEL ADAPTER
(USABLE WITH MANUAL TORQUE WRENCH)



Ref. 760CE

- 1 OT EQUATOR SCREWDRIVER/STEM FOR TORQUE CONTROLLER
CONTRAANGLE AND FOR MANUAL TORQUE WRENCH ADAPTER
READY FOR HOLDER (SQUARE 1,25mm)



Ref. 760CRD-US

- 1 MANUAL TORQUE WRENCH FOR OT Equator
TORQUE RANGE : 15-35 NCM - MAX . ALLOWED 50NCM



Ref. 124ICP

- 1 UNIVERSAL HANDLE



Ref.490IS

- 1 INSERTER/EXTRACTOR FOR SEEGER FOR FIXED PROSTHESIS
WITH TESTER, HANDLE AND CASE



Ref.190IS

- 1 INSERTER AND TESTER FOR FIXED PROSTHESIS SEEGER

MEDICAL STEEL T-BAR



Ref. G80A

- 4 MEDICAL STEEL JOINTS



Ref. T20A

- 4 MEDICAL STEEL CYLINDERS



Ref. T20Ai

- 4 MEDICAL STEEL RETENTIVE CYLINDERS

TITANIUM T-BAR



Ref. G80T

- 4 TITANIUM JOINTS



Ref. T20T

- 4 TITANIUM CYLINDERS



Ref. T20Ti

- 4 TITANIUM RETENTIVE CYLINDERS

TOOLS T-BAR



Ref. T-APLUS

- 2 T-BAR ADAPTERS

SPARE PARTS

	Ref. RMTSFEG	• 1 EXTRAGRADE TITANIUM ABUTMENT INCLINED AT 15° Ø 3,5mm
	Ref. RMTEG15SF	• 1 EXTRAGRADE TITANIUM ABUTMENT INCLINED AT 15° Ø 4mm
	Ref. RMTSEG	• 1 EXTRAGRADE TITANIUM ABUTMENT Ø 3,5mm
	Ref. RMT40SEG	• 1 EXTRAGRADE TITANIUM ABUTMENT Ø 4mm
	Ref. RMT50SEG	• 1 EXTRAGRADE TITANIUM ABUTMENT Ø 5mm
	Ref. RMMSEG	• 1 MINI EXTRAGRADE TITANIUM ABUTMENT Ø 3,5mm
	Ref. RMM40SEG	• 1 MINI EXTRAGRADE TITANIUM ABUTMENT Ø 4mm
	Ref. GCMR	• 1 CASTABLE SLEEVE FOR INCLINED TITANIUM ABUTMENT WITHOUT HOLE Ø 3,5mm
	Ref. GCMR4	• 1 CASTABLE SLEEVE FOR INCLINED TITANIUM ABUTMENT WITHOUT HOLE Ø 4mm
	Ref. GC	• 1 CASTABLE SLEEVE FOR TITANIUM ABUTMENT Ø 3,5mm
	Ref. GC4	• 1 CASTABLE SLEEVE FOR TITANIUM ABUTMENT Ø 4mm
	Ref. GC5	• 1 CASTABLE SLEEVE FOR TITANIUM ABUTMENT Ø 5mm
	Ref. GMM	• 1 CASTABLE SLEEVE FOR MINI TITANIUM ABUTMENT Ø 3,5mm
	Ref. GMM4	• 1 CASTABLE SLEEVE FOR MINI TITANIUM ABUTMENT Ø 4mm
	Ref. MTFP	• 1 TITANIUM ABUTMENT WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE WITH WIRE-REINFORCEMENT h 1,8 mm (ONLY AVAILABLE WITH Ø 3,5mm)
	Ref. MTB30FP	• 1 TITANIUM ABUTMENT WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE WITH WIRE-REINFORCEMENT h 3 mm (ONLY AVAILABLE WITH Ø 3,5mm)
	Ref. HR800	• 1 TITANIUM WIRE WITHOUT ELASTIC MEMORY Ø 1,0mm
	Ref. MEG40C	• 1 EXTRAGRADE CASTABLE ABUTMENT Ø 4mm
	Ref. RMCEG15SF	• 1 EXTRAGRADE CASTABLE ABUTMENT WITHOUT HOLE (ONLY AVAILABLE WITH Ø 4,0mm)
	Ref. SM6	• 6 WHITE SEEGER WITH HANDLE (STANDARD RETENTION)
	Ref. SRM6	• 6 PINK SEEGER WITH HANDLE (SOFT RETENTION)
	Ref. 022RGO	• 3 CASTABLE CONNECTION BAR
	Ref. FSS	• 1 BUR FOR SEEGER HOUSING POLISHING
	Ref. OC	• 1 OT CEM (COMPOSITE CEMENT: 1 Two-component syringe, 10 AUTOMIX TIPS)
	Ref. TOC	• 10 AUTOMIX TIPS

OT EQUATOR FOR IMPLANTS



Ref. 130

OT EQUATOR FOR IMPLANTS SET

SET CONTENT

- 1 OT EQUATOR TITANIUM + TiN ABUTMENT
- 1 STAINLESS STEEL HOUSING
- 1 PROTECTIVE DISK
- 4 RETENTIVE CAPS
(1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)



Ref. 131

OT EQUATOR AND SMARTBOX FOR IMPLANTS SET

SET CONTENT

- 1 OT EQUATOR TITANIUM + TiN ABUTMENT
- 1 ANODIZED SMART BOX HOUSING WITH BLACK POSITIONING CAP
- 1 PINK PROTECTIVE DISK
- 4 RETENTIVE CAPS
(1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)



Ref. 030

- 1 OT EQUATOR TITANIUM + TiN ABUTMENT

THE OT EQUATOR ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS. UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN ORDERING INDICATE: IMPLANT, DIAMETER AND CUFF HEIGHT. CUFF AVAILABLE TO 7mm.

OT EQUATOR FOR BARS



Ref. 160EQB

OT EQUATOR THREADED SLEEVE SYSTEM SET FOR BONDING

SET CONTENT

- 2 EQUATOR THREADED for titanium sleeve (1.6 mm THREAD)
- 2 THREADED SLEEVES (1.6 mm THREAD)
- 2 STAINLESS STEEL HOUSINGS
- 2 OT EQUATOR SPACERS
- 8 RETENTIVE CAPS
(2 EXTRA-SOFT, 2 SOFT, 2 STANDARD, 2 PROCESSING)

ELASTIC SEEGER PASSIVE BAR CONNECTIONS SYSTEM



Ref. 158ESA

ELASTIC SEEGER SET

SET CONTENT

- 1 CASTABLE CYLINDER FOR SEEGER
- 1 SELF-EXTRACTING SEEGER
- 1 TITANIUM SCREW FOR SELF-EXTRACTING SEEGER

CAPS ASSORTMENT SET



Ref. 192ECE

CAPS ASSORTMENT SET

SET CONTENT

- 1 STAINLESS STEEL HOUSING
- 1 PROCESSING CAP
- 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)
- 1 PROTECTIVE DISK

SMART BOX SET SELF-ALIGNING OT EQUATOR HOUSING



Ref. 335SBC

SMART BOX SET

SET CONTENT

- 1 SMART BOX HOUSING WITH BLACK POSITIONING CAP
- 1 PINK PROTECTIVE DISK
- 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)

OT EQUATOR CASTABLE
























Ref. 092ECQ

OT EQUATOR CASTABLE SET

SET CONTENT

- 2 SINGLE MALES
- 2 STAINLESS STEEL HOUSINGS
- 4 RETENTIVE CAPS (2 SOFT, 2 STANDARD)

SPARE PARTS OT EQUATOR

	Ref. 140CEV	• 4 RETENTIVE CAPS - VIOLET "STRONG" 2700g
	Ref. 140CET	• 4 RETENTIVE CAPS - CLEAR "STANDARD" 1800g
	Ref. 140CER	• 4 RETENTIVE CAPS - PINK "SOFT" 1200g
	Ref. 140CEG	• 4 RETENTIVE CAPS - YELLOW (EXTRA-SOFT 600g)
	Ref. 140CEN	• 4 PROCESSING CAPS - BLACK (FOR LABORATORY USE ONLY)
	Ref. 140CEDR8	• 4 RETENTIVE CAPS - AQUA CAPS <i>UNDERSIZED INTERNAL DIAMETER 2,2mm</i> (STRONG 1300g)
	Ref. 140CEADR8	• 4 RETENTIVE CAPS - ORANGE <i>UNDERSIZED INTERNAL DIAMETER 2,2mm</i> (VERY ELASTIC 350g)
	Ref. 141CAE	• 2 STAINLESS STEEL HOUSINGS
	Ref. 141CTE	• 2 TITANIUM HOUSINGS FUCHSIA ANODIZED
	Ref. 330SBE	• 1 SMART BOX HOUSING (WITH BLACK POSITIONING CAP)
	Ref. 335CSB	• 4 BLACK SMARTBOX (POSITIONING CAPS)
	Ref. 039SFE2	• 1 OT EQUATOR THREADED 2mm (DIRECT INTO 2mm THREAD)
	Ref. 339SFE	• 1 EQUATOR THREADED 1.6mm (FOR TITANIUM SLEEVE)
	Ref. 239GSF	• 1 TITANIUM SLEEVE (1.6mm THREAD)
	Ref. 339SSE	• 1 OT EQUATOR SPACER
	Ref. 159VEA	• 1 OT EQUATOR SEEGER (FOR SELF-EXTRACTING PEEK RING OR TITANIUM HEALING RING)
	Ref. 158VAE	• 1 TITANIUM SCREW FOR SELF-EXTRACTING SEEGER
	Ref. 158SAE	• 1 SELF-EXTRACTING WHITE SEEGER (FOR BAR LOCKING)
	Ref. 158CCE	• 1 CASTABLE CYLINDER height 2,5mm (FOR SEEGER)
	Ref. 158CCAE	• 1 HIGH CASTABLE CYLINDER height 3,5mm (FOR SEEGER)
	Ref. 151SS	• 2 CASTABLE MALES OT EQUATOR

ACCESSORIES OT EQUATOR



Ref. 144TLE

- 1 OT EQUATOR TITANIUM TRANSFER + TITANIUM SCREW



Ref. 144TTE

- 1 OT EQUATOR TITANIUM TRANSFER + TITANIUM SCREW



Ref. 144MTE

- 2 IMPRESSION COPINGS (FOR INDIVIDUAL TRAY)



Ref. 044CAIN

- 2 IMPRESSION COPINGS NORMAL SIZE - OT EQUATOR
(WITH INTERCHANGEABLE CAP FOR PICK UP IMPRESSION)



Ref. 144AE

- 2 LABORATORY ANALOGS ø 4mm

Ref. AI502

- 2 LABORATORY ANALOGS ø 5mm



Ref. 144AVC4

- 1 LABORATORY ANALOG WITH TITANIUM SCREW FOR CAD/CAM ø4mm

Ref. 144AVC5

- 1 LABORATORY ANALOG WITH TITANIUM SCREW FOR CAD/CAM ø5mm



Ref. 161SCME

- 4 OT EQUATOR SPACERS FOR HOUSING

OT EQUATOR TOOLS



Ref. 74AC01

- 1 PARALLELOMETER KEY "NORMAL SIZE"



Ref. 774CHE

- 1 OT EQUATOR SQUARE SCREW DRIVER + HOLDER
(SQUARE 1.25MM) (USABLE WITH RATCHET TORQUE CONTROL DEVICE)



Ref. 774HC

- 1 OT EQUATOR INTERCHANGEABLE HOLDER



Ref. 760CE

- 1 SQUARE DRIVER CONNECTOR FOR TORQUE
CONTROLLER CONTRA ANGLE (SQUARE 1.25mm)



Ref. 485IC

- 1 MULTIUSE CAPS INSERTER/EXTRACTOR TOOL



Ref. 185SIS

- 1 METAL INSERTION TOOL FOR SEEGER USE TOGETHER WITH 491EC



Ref. 491EC

- 1 MULTIUSE EXTRACTOR TOOL FOR CAPS
(WITH HOLE STEM HOLDER)



Ref. 760CRD-US

- 1 RATCHET TORQUE CONTROL DEVICE (FOR OT EQUATOR, SPHERO BLOCK
AND FLEX 15/35NCM STRENGTH - MAX 50NCM TORQUE, SUGGESTED 25NCM)

SINGLE THREADED SPHERES WITH THREADED SLEEVE TO BE BONDED NORMAL - MICRO



Ref. 139KSFN

SET SINGLE THREADED SPHERE NORMAL SIZE + TiN WITH SLEEVE TO BE BONDED SET CONTENT

- 2 TITANIUM NORMAL SINGLE THREADED SPHERE (HEX 1,3mm THREAD 1,6mm)
- 2 TITANIUM BOND-IN THREADED SLEEVE
- 2 S/STEEL SPACER for the THREADED SLEEVE (NORMAL sphere)



Ref. 139KSFM

SET SINGLE THREADED SPHERE MICRO SIZE + TiN WITH SLEEVE TO BE BONDED SET CONTENT

- 2 TITANIUM MICRO SINGLE THREADED SPHERE (HEX 0,9mm THREAD 1,6mm)
- 2 TITANIUM BOND-IN THREADED SLEEVE
- 2 S/STEEL SPACER for the THREADED SLEEVE (MICRO sphere)

SPARE PARTS



Ref. 239SFN

- 1 SINGLE THREADED SPHERE NORMAL SIZE + TiN FOR THREADED SLEEVE HEX 1,3mm THREAD 1,6mm



Ref. 239SFM

- 1 SINGLE THREADED SPHERE MICRO SIZE + TiN FOR THREADED SLEEVE HEX 0,9mm THREAD 1,6mm



Ref. 239GSF

- 1 TITANIUM BOND-IN THREADED SLEEVE FOR NORMAL AND MICRO - 1,6mm



Ref. 239SSN

- 1 STAINLESS STEEL SPACER FOR THREADED SLEEVE (NORMAL sphere)



Ref. 239SSM

- 1 STAINLESS STEEL SPACER FOR THREADED SLEEVE (Micro sphere)



Ref. 039SFN2

- 1 SINGLE THREADED SPHERE NORMAL SIZE + TiN HEX 1,3mm THREAD 2mm



Ref. 039SFM2

- 1 SINGLE THREADED SPHERE MICRO SIZE + TiN HEX 0,9mm THREAD 2mm



Ref. 760CBR

- 1 HEX DRIVER FOR THREADED SPHERE NORMAL HEX 1,27mm



Ref. 772CSF

- 1 HEX DRIVER FOR THREADED SPHERE MICRO HEX 0,9mm

OT CAP IMPRESSION COPINGS



Ref. 044CAIN

- 2 IMPRESSION COPINGS NORMAL SIZE Ø2,5mm WITH INTERCHANGEABLE CAP



Ref. 044CAI22

- 2 IMPRESSION COPINGS Ø2,25mm SPHERES WITH INTERCHANGEABLE CAP



Ref. 044CAIM

- 2 IMPRESSION COPINGS MICRO SIZE Ø1,8mm WITH INTERCHANGEABLE CAP

IMPLANTOLOGY TITANIUM ATTACHMENTS FOR OVERDENTURE

SPHERO FLEX - SPHERO BLOCK

Ref. 109



SPHERO FLEX TiN COATED

SET CONTENT

- 1 TITANIUM ABUTMENT WITH SWIVEL SPHERE (SPHERE Ø 2,5mm)
- 3 CAPS (Yellow - Pink - Clear)
- 1 STAINLESS STEEL HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

Ref. 002



SPHERO BLOCK NORMAL TiN COATED

SET CONTENT

- 1 TITANIUM ABUTMENT (SPHERE Ø 2,5 mm)
- 3 CAPS (Yellow - Pink - Clear)
- 1 STAINLESS STEEL HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

Ref. 003




SPHERO BLOCK MICRO TiN COATED

SET CONTENT

- 1 TITANIUM ABUTMENT (SPHERE Ø 1,8mm)
- 3 CAPS (Yellow - Pink - Clear)
- 1 STAINLESS STEEL HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

THE SPHERO FLEX, SPHERO BLOCK ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS. UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN ORDERING, PLEASE ALWAYS INDICATE: IMPLANT BRAND, IMPLANT DIAMETER AND CUFF HEIGHT. FOR SPHERO FLEX ABUTMENT CUFF AVAILABLE FROM 0,5 TO 7mm.



CUFF HEIGHT MEASURER


THIS TOOL IS USED TO DETERMINE THE CUFF HEIGHT OF IMPLANTS WITH INTERNAL OR EXTERNAL HEXAGON

Ref. 008MBG



SET CONTENT

- 1 CUFF HEIGHT SLIDER GAUGE
- 1 CUFF HEIGHT FIXED ROD GAUGE
- 1 SILICON RING DISPENSER
- 20 SILICON RINGS



TOOLS SPHERO FLEX - SPHERO BLOCK



Ref. 771CEF

- 1 KEY FOR SPHERO FLEX AND SPHERO BLOCK (NORMAL - Micro size) HEX 2,3mm




Ref. 760CBM

- 1 HEX SCREWDRIVER FOR TORQUE CONTROLLER CONTRA ANGLE HEX 1,3mm



Ref. 760CRD-US

- 1 RATCHET TORQUE CONTROL DEVICE
For Sphero block - flex and Ot Equator
15/35Ncm Strenght - Max 50Ncm torque, suggested 25Ncm.



Ref. 485IC

- 1 CAPS INSERTER/EXTRACTOR OT EQUATOR - NORMAL - MICRO



Ref. 008MDS

- 1 SPHERE MEASURING GAUGE
(TO CHECK THE SPHERICAL ATTACHMENT DIAMETERS: 1,6-1,8-2,25-2,5mm)




Ref. 044CAIN

- 2 IMPRESSION COPINGS NORMAL SIZE
Ø2,5mm WITH INTERCHANGEABLE CAP



Ref. 044CAIM

- 2 IMPRESSION COPINGS MICRO SIZE
Ø1,8mm WITH INTERCHANGEABLE CAP



Ref. 100AD

- DIRECTIONAL RINGS FOR SPHERO FLEX AND BLOCK
- 3 COLOURED RINGS WITH INCLINATION OF 0° - 7° - 14°



Ref. 044ABN

SPHERO FLEX - BLOCK NORMAL ANALOG SET SET CONTENT

- 1 SPHERO BLOCK-FLEX NORMAL ANALOG
- 1 BLACK NORMAL CAP FOR LABORATORY USE
- 3 DIRECTIONAL RINGS (0° - 7° - 14°)



Ref. 044ABM

SPHERO BLOCK MICRO ANALOG SET SET CONTENT

- 1 SPHERO BLOCK-MICRO ANALOG
- 1 BLACK MICRO CAP FOR LABORATORY USE
- 3 DIRECTIONAL RINGS (0° - 7° - 14°)

ANTI-UNSCREWING SYSTEM WITH AN ELASTIC DOWEL APPLICABLE TO

SPHERO FLEX-BLOCK - OT REVERSE AND DIFFERENT SCREWS. THE ELASTIC DOWEL IS APPLICABLE ONLY UPON REQUEST TO ALL SCREWS WITH A NOT INFERIOR M. 1.8 to screws - Sphero Flex - Sphero Block - OT Reverse attachments, etc.

Extra charge for
the dowel application



ELASTIC DOWEL



ATTACHMENTS OT CAP & OT BOX - OT STRATEGY - OT BAR OT VERTICAL - OT TECNO - TITANIUM SINGLE SPHERE MASTER STARTER SET FOR LABORATORY



SET CONTENT

OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (NORMAL and micro size)
- 4 castable spheres (NORMAL and micro size)
- 2 titanium single spheres (NORMAL and micro size)
+ sphere holder + 2 sliding base
- 4 castable OT CAP bars (NORMAL and micro size)
- 2 complete castable OT CAP TECNO bars (NORMAL and micro size)
- 4 castable OT BOX CLASSIC bars (upper + lower part)
NORMAL and micro size
- 2 castable OT BOX SPECIAL bars (NORMAL and micro size)
- 6 castable OT BOX connectors
- 4 castable OT BOX MONO housings (NORMAL and micro size)
- 4 plastic positioner rings (NORMAL and micro size)
- 16 OT CAP retentive caps (NORMAL and micro size)
Clear- standard retention - PINK soft retention
YELLOW extra soft retention - GREEN very elastic retention
- 4 TITAN CAP titanium caps (NORMAL and micro size)
- 4 BLACK CAPS- only for lab use - (NORMAL and micro size)
- 8 S/STEEL housings
- 4 ANALOGUES for plaster models (NORMAL and micro size)
- 4 protective disks
- 4 aluminium spacers (NORMAL and micro size)

OT STRATEGY:

- 4 castable spheres (2 standard and 2 high)
- 2 castable STEADY
- 6 retentive caps: 2 clear - 2 pink - 2 yellow for housing
- 2 black caps for housing only for lab use
- 2 s/steel housings
- 2 plastic positioner rings
- 6 retentive caps for duplication technique:
CLEAR - standard retention - PINK soft retention
YELLOW extra soft retention
- 2 black caps for duplication technique only for lab use
- 2 ANALOGUES for plaster models

OT BAR MULTIUSE:

- 2 castable bar
- 1 gingival connector
- 2 positioner clip A
- 2 positioner clip B
- 4 castable housings
- 4 clips - 2 PINK - soft retention - 2 YELLOW medium retention

Ref. 005SKLUS

OT VERTICAL:

- 2 castable male
- 2 castable STEADY
- 4 clips - 2 WHITE- standard retention - 2 GREEN soft retention
- 2 castable parallelometer keys
- 2 ceramic bars

TOOLS:

- 1 Pincers • 2 Reamer tools for caps (NORMAL and micro size)
- 1 Reamer tools for OT STRATEGY caps
- 2 Parallelometer keys for OT CAP (NORMAL and micro size)
- 1 Parallelometer key for OT STRATEGY size
- 2 Parallelometer keys for OT CAP TECNO (NORMAL and micro size)
- 1 Parallelometer key for OT BAR MULTIUSE
- 1 Blue Plastic "multiuse" insertion tool
(for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- 1 OT STRATEGY insertion tool
- 1 OT BAR clips insertion tool
- 1 OT VERTICAL clips insertion tool

ATTACHMENTS OT CAP & OT BOX - OT STRATEGY- OT BAR OT VERTICAL - OT UNILATERAL - OT EQUATOR

BASIC LABORATORY SET



SET CONTENT

OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (NORMAL and micro size)
- 8 castable spheres NORMAL/micro
- 2 castable OT CAP bars (NORMAL and micro size)
- 2 castable OT BOX CLASSIC bars (upper + lower part) NORMAL and micro size
- 2 castable OT BOX SPECIAL bars (NORMAL and micro size)
- 6 castable OT BOX connectors
- 4 castable OT BOX MONO housings (NORMAL and micro size)
- 8 plastic positioner rings (NORMAL and micro size)
- 13 OT CAP assorted retentive caps (NORMAL and micro size)
- PINK soft retention - YELLOW extra soft retention
- 13 BLACK CAPS- only for lab use - (NORMAL and micro size)
- 4 S/STEEL housings (NORMAL and micro size)

OT STRATEGY:

- 4 castable spheres (2 standard and 2 high)
- 2 castable STEADY
- 2 pink retentive caps
- 2 s/steel housings
- 2 plastic positioner rings
- 3 retentive caps for duplication technique:
- PINK soft retention - YELLOW extra soft retention
- 4 black caps for duplication technique only for lab use

OT EQUATOR:

- 2 single males
- 2 stainless steel housings
- 4 retentive caps (2 soft, 2 standard)
- 2 processing black caps

OT UNILATERAL:

- 1 unilateral castable attachment
- 1 uni-box castable
- 1 plastic positioner micro size
- 2 ot cap micro size: 1 pink - 1 black processing cap
- 2 ot strategy (for duplication technique): 1 yellow - 1 black processing cap
- 1 connector

OT BAR MULTIUSE:

- 1 castable bar
- 1 gingival connector
- 4 positioner clip (A - B)
- 2 castable housings
- 4 clips - 2 PINK - soft retention - 2 YELLOW medium retention

Ref. 005SKLBUS

OT VERTICAL:

- 2 castable male
- 2 castable STEADY
- 4 clips - 2 WHITE- standard retention - 2 GREEN soft retention
- 2 castable parallelometer keys
- 2 ceramic bars

TOOLS:

- 1 Tweezer
- 1 Parallelometer keys for OT CAP NORMAL
- 1 Parallelometer keys for OT CAP micro
- 1 Parallelometer key for OT STRATEGY
- 1 Parallelometer key for OT BAR MULTIUSE
- 1 Blue Plastic "multiuse" insertion tool
- (for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- 1 OT STRATEGY insertion tool
- 1 OT BAR clips insertion tool
- 1 OT VERTICAL clips insertion tool

OT CAP AND OT BOX ATTACHMENTS



NORMAL
Ref. 092CAN

MICRO
Ref. 092CAM

OT CAP NORMAL/MICRO

SET CONTENT

- 4 CASTABLE SINGLE SPHERES
- 4 PINK RETENTIVE CAPS
- 4 STAINLESS STEEL HOUSINGS
- 4 PLASTIC POSITIONERS



NORMAL
Ref. 196BCN

MICRO
Ref. 197BCM

OT CAP NORMAL/MICRO "ECONOMIC"

SET CONTENT

- 1 CASTABLE PRESECTIONED BAR
- 1 CASTABLE PRESECTIONED BEVELED BAR
- 4 CLEAR RETENTIVE CAPS (Standard retention)



NORMAL
Ref. 099BSN

MICRO
Ref. 099BSM

OT CAP & OT BOX MONO NORMAL/MICRO

SET CONTENT

- 2 CASTABLE PRESECTIONED BARS (1 straight + 1 beveled)
- 4 CLEAR RETENTIVE CAPS
- 4 CASTABLE MONO OT BOX
- 4 PLASTIC POSITIONERS



NORMAL
Ref. 058BSN

MICRO
Ref. 058BSM

OT BOX SPECIAL NORMAL/MICRO + CONNECTORS

SET CONTENT

- 2 OT BOX SPECIAL BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS



NORMAL
Ref. 153BCN

MICRO
Ref. 153BCM

OT BOX CLASSIC NORMAL/MICRO + CONNECTORS

SET CONTENT

- 2 UPPER BARS
- 2 LOWER BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

OT BOX LARGE OVERSIZED CASTABLE HOUSING FOR CHAIRSIDE POSITIONING OF CAPS



NORMAL
Ref. 253BLN

MICRO
Ref. 253BLM

OT BOX LARGE NORMAL/MICRO + CONNECTORS

SET CONTENT

- 2 OT BOX LARGE BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

OT STRATEGY ATTACHMENTS



Ref. 098SSS

OT STRATEGY CAPS FOR DUPLICATION TECHNIQUE

SET CONTENT

- 4 CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)
- 2 CASTABLE STEADY
- 4 RETENTIVE CAPS FOR DUPLICATION TECHNIQUE



Ref. 098CAL

OT STRATEGY CAPS FOR STAINLESS STEEL HOUSING TECHNIQUE

SET CONTENT

- 4 CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)
- 4 RETENTIVE CAPS FOR HOUSING TECHNIQUE (2 *SOFT*, 2 *STANDARD*)
- 4 STAINLESS STEEL HOUSINGS
- 2 PLASTIC POSITIONERS
- 2 CASTABLE STEADY

OT BAR MULTIUSE ATTACHMENTS



Ref. 021OBM

OT BAR MULTIUSE + CONNECTORS

SET CONTENT

- 2 CASTABLE BARS
- 4 POSITIONER CLIP A
- 4 POSITIONER CLIP B
- 4 CASTABLE BOXES
- 4 RETENTIVE PINK CLIPS
- 4 RETENTIVE YELLOW CLIPS
- 4 RETENTIVE RED CLIPS
- 2 CONNECTORS
- 1 GINGIVAL CONNECTOR

OT VERTICAL ATTACHMENTS



Ref. 071OBV

OT VERTICAL

SET CONTENT

- 4 CASTABLE MALES
- 2 CASTABLE STEADY
- 4 RETENTIVE WHITE CLIPS
- 2 RETENTIVE GREEN CLIPS
- 4 CASTABLE PARALLELOMETER KEYS + PIN

OT UNILATERAL EXTRA CORONAL ATTACHMENTS WITH COMBINED FUNCTION



Ref. 098UNL

OT UNILATERAL + CONNECTOR

SET CONTENT

- 1 UNILATERAL CASTABLE ATTACHMENT
- 1 UNI-BOX CASTABLE
- 1 PLASTIC POSITIONER MICRO SIZE
- 1 PINK CAP MICRO SIZE
- 1 STRATEGY YELLOW CAP (for duplication technique)
- 1 BLACK CAP MICRO SIZE (for laboratory use only)
- 1 STRATEGY BLACK CAP (for laboratory use only)
- 1 CONNECTOR



Ref. 156BCU

- 1 UNILATERAL CASTABLE ATTACHMENT



Ref. 156CCU

- 1 UNI-BOX CASTABLE

OT CAP TECNO

NORMAL

Ref. 093CTN

MICRO

Ref. 093CTM



OT CAP TECNO NORMAL/MICRO

SET CONTENT

- 2 PRE-ANGULATED CASTABLE EXTENSIONS
- 2 TITANIUM SINGLE THREADED SPHERES
- 2 TITANIUM BOND-IN THREADED SLEEVES
- 2 TRANSPARENT CAPS NORMAL/MICRO
- 2 PINK CAPS NORMAL/MICRO
- 2 YELLOW CAPS NORMAL/MICRO



Ref. 193BCA

- 2 PRE-ANGULATED CASTABLE EXTENSIONS

OT LOCK LOCKING PIN TITANIUM AND CASTABLE



Ref. 880CLT

OT LOCK SET TITANIUM

SET CONTENT

- 1 LOCKING PIN SET
- 1 HOUSING SHAPER
- 1 CERAMIC PIN



Ref. 880CLR

OT LOCK SET TITANIUM ADJUSTABLE

SET CONTENT

- 1 LOCKING PIN SET
- 1 HOUSING SHAPER
- 1 CERAMIC PIN
- 9 CASTABLE SPACER RINGS

OVERDENTURE ATTACHMENTS

S.P.L. PIVOTS - COPING COVER PROMOTIONAL SET FOR CLINIC



Ref. 011SPL

SET CONTENT

- 24 CASTABLE ASSORTED PIVOTS MICRO SIZE
- 40 CASTABLE ASSORTED PIVOTS NORMAL SIZE
- 12 TITANIUM MICRO PIVOTS
- 12 TITANIUM NORMAL PIVOTS
- 16 NORMAL CAPS (Clear - standard retention - PINK soft retention)
- 16 MICRO CAPS (Clear - standard retention - PINK soft retention)
- 20 PROTECTIVE DISKS
- 2 REGULATING TOOLS (1 NORMAL and 1 micro size)
- 1 BLUE PLASTIC "MULTIUSE" INSERTION TOOL
(FOR THE INSERTION OF THE CAPS INSIDE THE METAL HOUSINGS)
- 1 OT CAP insertion tool
- 4 NORMAL STAINLESS STEEL HOUSING
- 4 MICRO STAINLESS STEEL HOUSING
- 12 DIRECTIONAL RINGS (grey 0° - green 7° - red 14°)
- 1 TWEEZER



Lengths available:
10 - 9 - 7mm

Ref. 333PEK7

Ref. 333PEK9

Ref. 333PEK10

OT EQUATOR TITANIUM TIN COATED PIVOTS SET

LOW PROFILE HEAD PIVOTS

SET CONTENT

- 2 TITANIUM OT EQUATOR PIVOTS h7,9 or 10mm (of the same length chosen)
- 2 OT EQUATOR ANODIZED TITANIUM HOUSINGS
- 2 YELLOW EXTRA SOFT CAPS
- 2 PINK SOFT CAPS
- 2 CLEAR STANDARD CAPS
- 2 TRANSPARENT PROTECTIVE DISKS



Lengths available:
10 - 9 - 7mm

Ref. 336PTE7

Ref. 336PTE9

Ref. 336PTE10

OT EQUATOR TITANIUM TIN COATED PIVOTS

LOW PROFILE HEAD PIVOTS

SET CONTENT

- 1 TITANIUM OT EQUATOR PIVOT h7,9 or 10mm (of the length chosen)




Ref. 192ECE

OT EQUATOR CAPS ASSORTMENT SET

SET CONTENT

- 1 STAINLESS STEEL HOUSING
- 1 PROCESSING CAP
- 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)
- 1 PROTECTIVE DISK



Ref. 033PSF7R

Ref. 033PSF9R

Ref. 033PSF10R

PIVOT FLEX TIN COATED

TITANIUM PIVOT WITH SWIVEL SPHERE - NORMAL size Adapted for COPING COVER

SET CONTENT

- 1 TITANIUM PIVOT WITH SWIVEL SPHERE
- 1 STAINLESS STEEL HOUSING
- 2 PINK CAPS NORMAL Size - Soft retention
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

Disponibile
nelle lunghezze:
10 - 9 - 7mm
Sfera NORMO ø2,5
Sfera MICRO ø1,8



Ref. 036PTN7

Ref. 036PTN9

Ref. 036PTN10

TITANIUM PIVOTS SET NORMAL SIZE

SET CONTENT

- 5 TITANIUM PIVOTS (Adapted for COPING COVER)

Ref. 036PTM7

Ref. 036PTM9

Ref. 036PTM10

TITANIUM PIVOTS SET MICRO SIZE

SET CONTENT


- 5 TITANIUM PIVOTS (Adapted for COPING COVER)

NORMAL

Ref. 192ACN

MICRO

Ref. 192ACM



CAPS ASSORTMENT SET OT CAP NORMAL/MICRO

SET CONTENT

- 1 STAINLESS STEEL HOUSING
- 1 PROCESSING CAP
- 4 RETENTIVE CAPS (1 GREEN, 1 YELLOW, 1 PINK, 1 CLEAR)
- 1 PROTECTIVE DISK



Ref. 034PRK

OT REVERSE 3 OVERDENTURE POSTS

FOR ROOTS, AVAILABLE IN ONE STANDARD SIZE

SET CONTENT

- 2 ROOT POSTS IN TITANIUM + TIN
- 2 RETENTIVE MALE ATTACHMENT TITANIUM + NYLON
- 2 PLASTIC MANUAL TOOLS
- 2 STAINLESS STEEL HOUSINGS
- 2 PINK PROTECTIVE DISKS

Ø del perno sottosfera:
2,2 2 1,8 1,7 1,4mm



Ref. 010PSP

CASTABLE PIVOTS FOR CANAL IMPRESSIONS

SET CONTENT

- 20 CASTABLE PIVOTS NORMAL SIZE

Lengths available:
14 - 12 - 10 - 9 - 7mm



Ref. 012PSM

CASTABLE PIVOTS FOR CANAL IMPRESSIONS

SET CONTENT

- 20 CASTABLE PIVOTS MICRO SIZE

Lengths available:
10 - 9 - 7mm



Ref. A01MOG

MOOSER BUR FOR PIVOTS OF 7 - 9 - 10mm



Ref. A03MOB

MOOSER BUR FOR PIVOTS OF 12 - 14mm



Ref. 121FD

OT REVERSE 3 DIAMOND SHAPED BUR FOR ROOTS



Ref. 100AD

DIRECTIONAL RINGS FOR SPHERO FLEX AND BLOCK

- 3 COLOURED RINGS WITH INCLINATION OF 0° - 7° - 14°



Ref. 485IC

- 1 EXTRACTOR/INSERTER TOOL FOR CAPS OT EQUATOR - NORMAL - MICRO

SPARE PARTS OT CAP SYSTEM NORMAL Ø 2,5 mm



Ref. 055SCN

• 4 CASTABLE SINGLE SPHERES NORMAL SIZE

Ref. 151BAN

• 2 CASTABLE PRESECTIONED BEVELED BARS NORMAL SIZE

Ref. 150BPN

• 2 CASTABLE PRESECTIONED BARS NORMAL SIZE



Ref. 048CAN

• 6 SILVER CAPS NORMAL SIZE (EXTRA-RESILIENT)
(ELASTIC AND GUMMY) - 350g



Ref. 048CON

• 6 GOLD CAPS NORMAL SIZE (EXTRA-RESILIENT)
(SLIGHTLY ELASTIC) - 500g



Ref. 049PCN

• 6 GREEN CAPS NORMAL SIZE
(VERY ELASTIC RETENTION) - 350g



Ref. 060CRN AY

• 6 YELLOW CAPS NORMAL SIZE
(EXTRA SOFT RETENTION) - 500g



Ref. 040CRN SN

• 6 PINK CAPS NORMAL SIZE
(SOFT RETENTION) - 900g



Ref. 040CRN

• 6 CLEAR CAPS NORMAL SIZE
(STANDARD RETENTION) - 1300g



Ref. 040TCN

• 2 TITAN CAPS NORMAL SIZE



Ref. 043CLN

• 6 BLACK CAPS NORMAL SIZE (FOR LABORATORY USE ONLY)



Ref. 049PCNDR8

• 6 ORANGE CAPS NORMAL SIZE (VERY ELASTIC RETENTION)
(UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 350g



Ref. 060CRNAYDR8

• 6 YELLOW CAPS NORMAL SIZE (EXTRA SOFT RETENTION)
(UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 500g



Ref. 040CRNSNDR8

• 6 PINK CAPS NORMAL SIZE (SOFT RETENTION)
(UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 900g



Ref. 040CRNDR8

• 6 AQUA CAPS NORMAL SIZE (STRONG RETENTION)
(UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 1300g



Ref. 042AMN

• 2 MONO OT BOX NORMAL SIZE
(CASTABLE HOUSINGS FOR CAPS)



Ref. 041CAN

• 2 STAINLESS STEEL HOUSINGS NORMAL SIZE



Ref. 041CTN

• 2 TITANIUM HOUSINGS NORMAL SIZE



Ref. 101PPN

• 2 PLASTIC POSITIONERS NORMAL SIZE



Ref. 044PPN

• 2 NORMAL SPHERE ANALOGS

SPARE PARTS OT CAP SYSTEM MICRO Ø 1,8 mm



Ref. 055SCM

• 4 CASTABLE SPHERES MICRO SIZE

Ref. 151BAM

• 2 CASTABLE PRESECTIONED BEVELED BARS MICRO SIZE

Ref. 150BPM

• 2 CASTABLE PRESECTIONED BARS MICRO SIZE

Ref. 048CAM

• 6 SILVER CAPS MICRO SIZE (EXTRA-RESILIENT)
(ELASTIC AND GUMMY) - 200g

Ref. 048COM

• 6 GOLD CAPS MICRO SIZE (EXTRA-RESILIENT)
(SLIGHTLY ELASTIC) - 450g

Ref. 049PCM

• 6 GREEN CAPS MICRO SIZE
(VERY ELASTIC RETENTION) - 200g

Ref. 060CRM AY

• 6 YELLOW CAPS MICRO SIZE
(EXTRA SOFT RETENTION) - 450g

Ref. 040CRM SN

• 6 PINK CAPS MICRO SIZE
(SOFT RETENTION) - 800g

Ref. 040CRM

• 6 CLEAR CAPS MICRO SIZE
(STANDARD RETENTION) - 1100g

Ref. 040TCM

• 2 TITAN CAPS MICRO SIZE

Ref. 043CLM

• 6 BLACK CAPS MICRO SIZE (FOR LABORATORY USE ONLY)

Ref. 049PCMDR8

• 6 ORANGE CAPS MICRO SIZE (VERY ELASTIC RETENTION)
(UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø1,6mm) - 200g

Ref. 040CRMDR8

• 6 AQUA CAPS MICRO SIZE (STRONG RETENTION)
(UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø1,6mm) - 1100g

Ref. 042AMM

• 2 MONO OT BOX MICRO SIZE
(CASTABLE HOUSINGS FOR CAPS)

Ref. 041CAM

• 2 STAINLESS STEEL HOUSINGS MICRO SIZE

Ref. 041CTM

• 2 TITANIUM HOUSINGS MICRO SIZE

Ref. 101PPM

• 2 PLASTIC POSITIONERS MICRO SIZE

Ref. 044PPM

• 2 MICRO SPHERE ANALOGS

OT CAP ASSORTMENT SET



KIT NORMAL

Ref. 074KAN

KIT MICRO

Ref. 074KAM

SET CONTENT

- 8 CLEAR CAPS - STANDARD RETENTION
- 8 PINK CAPS - SOFT RETENTION
- 6 YELLOW CAPS - EXTRA SOFT RETENTION
- 6 GREEN CAPS - VERY ELASTIC RETENTION
- 6 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC
- 6 SILVER CAPS - (EXTRA RESILIENT) ELASTIC AND GUMMY
- 4 LIGHT GREEN CAPS - STANDARD (FOR WORN OUT SPHERE)
- 4 ORANGE CAPS - ELASTIC (FOR WORN OUT SPHERE)
- 4 BLACK CAPS - LABORATORY USE ONLY
- 2 TITAN CAPS



Ref. 074KNM

ASSORTMENT SET NORMAL & MICRO

SET CONTAINS

- 4+4 CLEAR CAPS - STANDARD RETENTION
- 4+4 PINK CAPS - SOFT RETENTION
- 2+2 YELLOW CAPS - EXTRA SOFT RETENTION
- 3+3 GREEN CAPS - VERY ELASTIC RETENTION
- 3+3 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC
- 3+3 SILVER CAPS - (EXTRA RESILIENT) ELASTIC AND GUMMY
- 2+2 LIGHT GREEN CAPS - STANDARD (FOR WORN OUT SPHERE)
- 2+2 ORANGE CAPS - ELASTIC (FOR WORN OUT SPHERE)
- 2+2 BLACK CAPS - LABORATORY USE ONLY
- 1+1 TITAN CAPS

OT CAP ASSORTED RETENTIVE CAPS



NORMAL

Ref. 064ACN

MICRO

Ref. 064ACM

SET CONTENT

- 6 GREEN CAPS - VERY ELASTIC RETENTION
- 6 YELLOW CAPS - EXTRA SOFT RETENTION
- 6 PINK CAPS - SOFT RETENTION
- 6 CLEAR CAPS - STANDARD RETENTION

OT CAP SYSTEM ATTACHMENTS



NORMAL

Ref. 038STN

TITANIUM SINGLE SPHERE NORMAL SIZE + TiN

FOR LASER SOLDERING - FOR BONDING

SET CONTENT

- 2 SINGLE TITANIUM SPHERES
- 2 PINK CAPS - SOFT RETENTION
- 1 SPHERE HOLDER
- 2 CASTABLE SLIDING BASES



MICRO

Ref. 038STM











TITANIUM SINGLE SPHERE MICRO SIZE + TiN

FOR LASER SOLDERING - FOR BONDING

SET CONTENT

- 2 SINGLE TITANIUM SPHERES
- 2 PINK CAPS - SOFT RETENTION
- 1 SPHERE HOLDER
- 2 CASTABLE SLIDING BASES

OT CAP NORMAL-MICRO TOOLS

	Ref. 485IC	• 1 EXTRACTOR/INSERTER TOOL FOR CAPS OT EQUATOR - OT CAP NORMAL/MICRO
	Ref. 74AC01	• 1 PARALLELOMETER MANDREL NORMAL/OT EQUATOR
	Ref. 74AC02	• 1 PARALLELOMETER MANDREL MICRO
	Ref. 74AC03	• 1 PARALLELOMETER MANDREL TECNO NORMAL AND MICRO SIZE (ONLY FOR OT CAP TECNO)
	Ref. 491EC	• 1 MULTIUSE REMOVAL TOOL FOR RETENTIVE CAPS
	Ref. 008MDS	• 1 SPHERE MEASURING GAUGE (TO CHECK THE SPHERICAL ATTACHMENT DIAMETERS: 1,6-1,8-2,25-2,5mm)
	Ref. 080RCN	• 1 REAMER TOOL FOR OT CAPS - NORMAL SIZE
	Ref. 080RCM	• 1 REAMER TOOL FOR OT CAPS - MICRO SIZE
	Ref. 082ATN	• 1 TOOL FOR TESTING CAP RETENTION - NORMAL SIZE
	Ref. 083ATM	• 1 TOOL FOR TESTING CAP RETENTION - MICRO SIZE

NORMAL - MICRO ACCESSORIES AND Ø2.2mm

	Ref. 155CC	• 6 CONNECTORS
	Ref. 187CRS	• 2 CONCAVE SPHERES AVAILABLE: Ø1.8 - 2.2 - 2.5mm
	Ref. 187CRE	• 2 CONCAVE OT EQUATOR
	Ref. 189SRS	• 2 SOLID SPHERES AVAILABLE: Ø1.8mm
	Ref. 100PD	• 10 PROTECTIVE DISKS (internal diameter 1,6mm)
	Ref. 100PDR	• 10 PINK PROTECTIVE DISKS (internal diameter 2,9mm)
	Ref. 161SCMN	• 4 OT CAP NORMAL/MICRO SPACERS FOR HOUSING
	Ref. 161SCMM	
	Ref. 038GCS	• 2 CASTABLE SLIDING BASES
	Ref. 044PP22	• 2 MICRO ANALOGS Ø2.2mm

OT STRATEGY SPARE PARTS



Ref. 154PCS

- 4 OT STRATEGY STANDARD MALES



Ref. 054PLS

- 4 OT STRATEGY HIGH MALES



Ref. 154STY

- 2 CASTABLE STEADY



Ref. 047CSD AY

- 4 STRATEGY YELLOW CAPS FOR DUPLICATION TECHNIQUE
(EXTRA SOFT RETENTION) - 500g



Ref. 047CSD SN

- 4 STRATEGY PINK CAPS FOR DUPLICATION TECHNIQUE
(SOFT RETENTION) - 900g



Ref. 047CSD

- 4 STRATEGY CLEAR CAPS FOR DUPLICATION TECHNIQUE
(STANDARD RETENTION) - 1300g



Ref. 047CLS

- 4 STRATEGY BLACK CAPS FOR DUPLICATION TECHNIQUE
(FOR LABORATORY USE ONLY)



Ref. 047CSDDR8

- 4 AQUA STRATEGY CAPS FOR DUPLICATION TECHNIQUE
(UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE)



Ref. 045CSM AY

- 4 STRATEGY YELLOW CAPS FOR S/STEEL HOUSING
(EXTRA SOFT RETENTION) - 500g
(TO BE USED ONLY WITH HOUSINGS)



Ref. 045CSM SN

- 4 STRATEGY PINK CAPS FOR S/STEEL HOUSING
(SOFT RETENTION) - 900g
(TO BE USED ONLY WITH HOUSINGS)



Ref. 045CSM

- 4 STRATEGY CLEAR CAPS FOR S/STEEL HOUSING
(STANDARD RETENTION) - 1300g
(TO BE USED ONLY WITH HOUSINGS)



Ref. 045CLS

- 4 STRATEGY BLACK CAPS FOR S/STEEL HOUSING
(FOR LABORATORY USE ONLY)



Ref. 045CSMDR8

- 4 AQUA STRATEGY CAPS FOR S/STEEL HOUSING
(UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE)



Ref. 046CSM

- 2 STRATEGY STAINLESS STEEL HOUSINGS



Ref. 101OTS

- 2 STRATEGY PLASTIC POSITIONERS
(FOR SOLDERING STAINLESS STEEL HOUSINGS)

OT STRATEGY ASSORTMENT CAP SET FOR DUPLICATION TECHNIQUE

SET CONTENT

- 4 YELLOW CAPS - EXTRA SOFT RETENTION
- 4 PINK CAPS - SOFT RETENTION
- 4 CLEAR CAPS - STANDARD RETENTION

Ref. 047ACS

OT STRATEGY ASSORTMENT CAP SET FOR S/STEEL HOUSINGS

SET CONTENT

- 4 YELLOW CAPS - EXTRA SOFT RETENTION
- 4 PINK CAPS - SOFT RETENTION
- 4 CLEAR CAPS - STANDARD RETENTION

Ref. 045ACS

Ref. 044PPS

- 2 STRATEGY ANALOGS

Ref. 486ICS

- 1 INSERTER/EXTRACTOR TOOL FOR STRATEGY CAPS

Ref. 75AC04

- 1 PARALLELOMETER MANDREL FOR STRATEGY CAPS

Ref. 081RCS

- 1 REAMER TOOL FOR STRATEGY CAPS



OT BAR MULTIUSE SPARE PARTS

Ref. 022OBB

- 2 CASTABLE BARS OT BAR MULTIUSE

Ref. 025CPB

- 4 CASTABLE BOXES OT BAR MULTIUSE

Ref. 023CPA

- 4 POSITIONER CLIPS A

Ref. 024CPB

- 4 POSITIONER CLIPS B

Ref. 027CRG

- 4 YELLOW CLIPS (MEDIUM RETENTION 1,8Kg)

Ref. 026CRR

- 4 PINK CLIPS (SOFT RETENTION 1,2Kg)

Ref. 027CRS

- 4 RED CLIPS (EXTRA SOFT RETENTION 0,6Kg)

Ref. 022RGO


- 3 CASTABLE CONNECTORS

Ref. 429IOBM

- 1 INSERTER/EXTRACTOR TOOL FOR OT BAR MULTIUSE CLIPS

Ref. 028OCP

- 1 PARALLELOMETER KEY FOR OT BAR MULTIUSE



OT VERTICAL SPARE PARTS



Ref. 072MCV

- 2 CASTABLE MALES
- 2 PIN KEYS

Ref. 072SCV

- 2 CASTABLE STEADY

Ref. 072CBV

- 4 WHITE CLIPS (STANDARD RETENTION)

Ref. 072CVV

- 4 GREEN CLIPS (SOFT RETENTION)

Ref. 472ICV

- 1 INSERTION TOOL FOR OT VERTICAL CLIPS

OT REVERSE 3 SPARE PARTS



Ref. 120MR2

- 2 RETENTIVE MALES IN TITANIUM + NYLON (COLOR CLEAR/ STANDARD RETENTION)

Ref. 041CAM

- 2 S/STEEL HOUSINGS MICRO SIZE

Ref. 121AAR

- 2 INOX ROOTS ANALOGS

Ref. 121AM

- 2 PLASTIC MANUAL TOOLS

Ref. 121PR

- 2 PLASTIC POSITIONERS FOR DUPLICATION

Ref. 100PDR

- 10 PINK PROTECTIVE DISKS

Ref. 124IMR

- 1 METAL INSERTING TOOL FOR OT REVERSE/3 MALE

Ref. 72CBR

- 1 S/STEEL PARALLELOMETER KEY FOR CASTABLE BAR

Ref. 72CFR

- 1 BRASS PARALLELOMETER KEY FOR TITANIUM FEMALE

OT LOCK SPARE PARTS



Ref. 882CS

- 1 SLIDING KEY
(ONLY FOR TITANIUM - CASTABLE LOCKING PIN)



Ref. 882AF

- 1 FLEXIBLE RING



Ref. 882CB

- 1 LOCKING KEY



Ref. 882CPT

- 1 TITANIUM PREFABRICATED HOUSING
(ONLY FOR TITANIUM - TITANIUM ADJUSTABLE LOCKING PIN)



Ref. 882PPC

- 1 CERAMIC PIN



Ref. 882CLSCT

- 1 ADJUSTABLE SLIDING KEY
(ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)



Ref. 882DCC

- 9 CASTABLE SPACER RINGS
(ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)

OT LOCK ACCESSORIES



Ref. 882CG

- 1 CONICAL GUIDE



Ref. 882CAS

- 1 UNLOCKING TOOL

SPHERE REPAIR SETS



Ø AVAILABLE:
1,8 - 2,25 - 2,5

Ref. 087CRS

CONCAVE REPAIR SPHERE

SET CONTENT

- 2 CONCAVE SPHERES IN TITANIUM - TiN COATED
- 2 PINK CAPS SOFT RETENTION
- 1 INSERTING TOOL
- 1 GAUGE AND STRIP HOLDER



Ref. 087CRE

CONCAVE REPAIR OT EQUATOR

SET CONTENT

- 2 CONCAVE OT EQUATOR IN TITANIUM - TiN COATED
- 2 PINK CAPS SOFT RETENTION
- 1 INSERTING TOOL
- 1 GAUGE AND STRIP HOLDER



Ref. 089SRS

SOLID REPAIR SPHERE

SET CONTENT

- 2 SOLID SPHERES 1.8mm IN TITANIUM - TiN COATED
- 2 PINK CAPS SOFT RETENTION
- 2 PROTECTIVE DISKS
- 1 KEY

DENTURE BOXES



Ref. 0P6

SET CONTENT


- 6 ASSORTED BOXES, COLOURS: BLUE - LIGHT BLUE - WHITE



Ref. 100AS

- 30 SILICONE RINGS

SMALL PARALLELOMETER



Ref. 00PB

- 1 **SMALL PARALLELOMETER DEVICE**
WITH MODEL HOLDER BASE
(FOR COURSES, ETC.)


OT LEVER

UNIVERSAL LEVER



Ref. 761OTL

- 2 **OT LEVER**
FOR REMOVAL OF PROSTHESES AND ALIGNERS



Ref. BOL

OT LEVER BOX

SET CONTENT

- 6 OT LEVER PACKAGES
- 1 OT LEVER PACKAGING DISPLAY STAND



OT CEM

COMPOSITE CEMENT



Ref. OC

SET CONTENT

- 1 TWO-COMPONENT SYRINGE 5ml
- 10 AUTOMIX TIPS




Ref. TOC

- 10 AUTOMIX TIPS

BROKEN SCREW EXTRACTOR KIT

BROKEN SCREW EXTRACTOR SETS ARE MANUFACTURED FOR MOST IMPLANT BRANDS



Ref. 680

SET CONTENT

- 1 CLAW REAMER BUR
- 1 REVERSE CUTTING BUR
- 1 MANUAL EXTRACTOR
- 1 CENTERING DEVICE

Ref. 680FS

- 1 REVERSE CUTTING BUR

Ref. 680FA*

- 1 CLAW REAMER BUR

*Specify type of implant and diameter

THE BROKEN SCREW EXTRACTOR SET MAY HAVE SPACE PROBLEMS IN THE REAR SECTION, WHICH COULD AFFECT ITS PERFORMANCE. IT'S IMPORTANT TO ANALYZE THE SITUATION BEFORE USING THE SET.

Scientific training



Different approaches to dentistry are proposed, from traditional techniques to the most advanced digital solutions.

Live courses




The surgical-prosthetic process in all its phases, an excellent tool to facilitate learning new protocols.

RHEIN83

All on
40
1983-2023

Research and Development




Innovation and research is developed with the collaboration of numerous universities, both Italian and foreign.


Academy Rhein83




Conferences, events and practical courses, welcoming professionals from all over the world for training in fixed and removable prosthetics.



**View Rhein83
calendar courses**



**Consult the Rhein83's
scientific literature**



**Watch Rhein83's
clinic videos**



DIRECT SERVICES

A telephone service is available during office hours,
it furnishes direct answers for technical questions.
You can address to Rhein83 for knowing the address of foreign distributors

Tel. +39 051 244510 - +39 051 244396 - Fax +39 051 245238

Rhein83 website
www.rhein83.com

General informations
nick@rhein83.it

Sales department
international.orders@rhein83.it

Technical support
direzione.tecnica@rhein83.it

Marketing department
direzione@rhein83.it



RHEIN83

Via ZAGO, 10/ABC
40128 - BOLOGNA (ITALY)

Tel. (+39) 051 244510 - (+39) 051 244396
Fax (+39) 051 245238

MOD. D395SP REV.01 del 12/09/2023

